

# MICROBIAL SEALANT



## MICROBIAL SEALANT

### HALYARD\* INTEGUSEAL\* Microbial Sealant

- Ergonomic Single-Use Applicator
- Evenly Covers Incision Area In A Single Stroke
- Widely-Used Cyanoacrylate Forms A Safe, Fast-Drying, Breathable Barrier
- Compatible With Most Common Surgical Techniques And Materials
- Need Not Be Removed For Skin Closure
- Can Be Used With A Variety Of Preps And Incise Drapes
- Nonflammable, Latex-Free, Non-Irritating And Non-Cytotoxic

Code	Description	Coverage Area (approx)	Example of Procedures	Packaging
33736	IS 100	10" X 10" 25cm X 25cm	Cranial and Spinal Procedures, Laparotomies, Hernia Repair, Thyroid Surgery, Cardiac Catheterisation	Applicators/Box = 10 Boxes/Case = 2 Applicators/Case = 20
33737	IS 200	20" X 10" 50cm X 25cm	Cholecystectomy, Open Heart, Lung Surgery, Total Joint Replacement	Applicators/Box = 10 Boxes/Case = 2 Applicators/Case = 20