

RESPIRATORY HEALTH

CLOSED SUCTION SYSTEMS / ACCESSORIES

HALYARD* MULTI-ACCESS PORT CLOSED SUCTION SYSTEM

- Rotating Manifold
- Single Use Cassette
- Turbulent Cleaning Chamber
- Sleeve Tether
- MDI Drug Delivery
- Replacement Catheters

Code	Description	Packaging
240128-5	HALYARD* Multi-Access Port Closed Suction System for Adults, Endotracheal, MDI, 12 F	20/Case
240148-5	HALYARD* Multi-Access Port Closed Suction System for Adults, Endotracheal, MDI, 14 F	20/Case
2401238-5	HALYARD* Multi-Access Port Closed Suction System for Adults, Tracheostomy, DSY, MDI, 12 F	20/Case
2401438-5	HALYARD* Multi-Access Port Closed Suction System for Adults, Tracheostomy, DSY, MDI, 14 F	20/Case
240125-5	HALYARD* Multi-Access Port Replacement Catheter for Adults, Endotracheal, 12 F	20/Case
240145-5	HALYARD* Multi-Access Port Replacement Catheter for Adults, Endotracheal, 14 F	20/Case
2401253-5	HALYARD* Multi-Access Port Replacement Catheter for Adults, Tracheostomy, 12 F	20/Case
2401453-5	HALYARD* Multi-Access Port Replacement Catheter for Adults, Tracheostomy, 14 F	20/Case
24013-5	HALYARD* Seal Cassette	20/Case

NOTE: Please see package inserts for directions for use for each of these products.


RESPIRATORY HEALTH

TURBO ADULT CLOSED SUCTION SYSTEMS

HALYARD* Adult TRACH CARE* 72 Hour Double Swivel Elbows


- Gamma Sterilised
- One-Way Lavage / Irrigation Port
- Dual Swivels Provides Flexibility Without Connectors
- Tactile Sleeve for Better Feel

- Numbered Catheter Allowing the Clinician to Limit the Distance of Catheter Advance
- Tracheostomy Sizes



Code	Outer Diameter	Length	Description	Packaging
227-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow Endotracheal	20/Case
22703-5	14 Fr / 4.6mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case
22716-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow Endotracheal	20/Case
2271603-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case

HALYARD* Adult TRACH CARE* 72 Hour T-Pieces



Code	Outer Diameter	Length	Description	Packaging
22705-5	14 Fr / 4.6mm	21.3" / 54cm	T-Piece Endotracheal	20/Case
2270135-5	14 Fr / 4.6mm	12" / 30.5cm	T-Piece Tracheostomy	20/Case
22715-5	12 Fr / 4mm	21.3" / 54cm	T-Piece Endotracheal	20/Case
227135-5	12 Fr / 4mm	12" / 30.5cm	T-Piece Tracheostomy	20/Case

NOTE: Please see package inserts for directions for use for each of these products.


RESPIRATORY HEALTH

CLOSED SUCTION SYSTEMS / ACCESSORIES

KIMVENT* CLOSED SUCTION SYSTEMS

HALYARD* Adult TRACH CARE* 24 Hour Double Swivel Elbows


- Gamma Sterilised
- One-Way Lavage / Irrigation Port
- Dual Swivels Provides Flexibility Without Connectors
- Tactile Sleeve for Better Feel
- Numbered Catheter Allowing the Clinician to Limit the Distance of Catheter Advance
- Tracheostomy Sizes



Code	Outer Diameter	Length	Description	Packaging
2160-5	12 Fr / 4mm	21.3" / 54cm	Double Swivel Elbow	20/Case
21603-5	12 Fr / 4mm	12" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case
2210-5	14 Fr / 4.6mm	21.3" / 54cm	Double Swivel Elbow	20/Case
22103-5	14 Fr / 4.6mm	21.3" / 30.5cm	Double Swivel Elbow Tracheostomy	20/Case
2260-5	16 Fr / 5.3mm	21.3" / 54cm	Double Swivel Elbow	20/Case

HALYARD* Adult TRACH CARE* 24 Hour T-Pieces

- Gamma Sterilised
- Lavage / Irrigation Port
- Locking Thumb Port
- T-Piece
- Multiple French Sizes



Code	Outer Diameter	Length	Description	Packaging
2105-5	10 Fr / 3.3mm	21.3" / 54cm	Clear T-Piece	20/Case
2155-5	12 Fr / 4mm	21.3" / 54cm	Clear T-Piece	20/Case
2205-5	14 Fr / 4.6mm	21.3" / 54cm	Clear T-Piece	20/Case
220135-5	14 Fr / 4.6mm	12" / 30.5cm	Clear T-Piece Tracheostomy Swivel Adapter	20/Case
2255-5	16 Fr / 5.3mm	21.3" / 54cm	Clear T-Piece	20/Case
25-2210-5	14 Fr / 4.6mm	21.3" / 54cm	24hr Double Swivel Elbow With Safety Lavage Port	20/Case

NOTE: Please see package inserts for directions for use for each of these products.

RESPIRATORY HEALTH

PAEDIATRIC / NEONATAL CLOSED SUCTION SYSTEMS

HALYARD* Neonatal and Paediatric TRACH CARE* "Y"



Code	Outer Diameter	Length	Description	Packaging
195-5	5 Fr / 1.66mm	12" / 30.5cm	"Y" Endotracheal Tube Adapters (2.0mm, 2.5mm)	20/Case
196-5	6 Fr / 2mm	12" / 30.5cm	"Y" Endotracheal Tube Adapters (2.5mm, 3.0mm, 3.5mm)	20/Case
197-5	7 Fr / 2.3mm	12" / 30.5cm	"Y" Endotracheal Tube Adapters (3.0mm, 3.5mm, 4mm)	20/Case
198-5	8 Fr / 2.6mm	12" / 30.5cm	"Y" Endotracheal Tube Adapters (3.5mm, 4.0mm, 4.5mm)	20/Case
1910-5	10 Fr / 3.3mm	16" / 40.5cm	"Y" Endotracheal Tube Adapters (4.5mm, 5.0mm, 5.5mm)	20/Case
1912-5	12 Fr / 4mm	16" / 40.5cm	"Y" Endotracheal Tube Adapters (5.5mm, 6.0mm, 6.5mm)	20/Case

HALYARD* Neonatal and Paediatric MAC Catheters



Code	Outer Diameter	"Y" Size / Plug Size	Description	Packaging
1900	5 Fr / 1.66mm	2.5mm / 7.5mm, 3.0mm / 7.5mm, 3.5mm / 7.5mm	MAC Catheter	5/Case
1940	5 Fr / 1.66mm	4.0mm / 7.5mm	MAC Catheter	5/Case

NOTE: Please see package inserts for directions for use for each of these products.