


# RESPIRATORY HEALTH

## **HALYARD\* MICROCUFF\*** **ADULT & PAEDIATRIC ENDOTRACHEAL TUBES**


### **HALYARD\* MICROCUFF\* Endotracheal Tube, Oral/Nasal Magill**

- Ultra-low Cuff Pressure Can Reduce the Risk of Tracheal Trauma
- Short, Cylindrical Cuff Placed Near the Tracheal Tube Tip
- Four Precision Bands to Facilitate and Confirm Optimal Tube Placement
- No Murphy Eye
- Oral/Nasal Magill and Pre-Curved Shapes



Code	Inner Diameter	Description	Packaging
35111	3.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35112	3.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35113	4.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35114	4.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35115	5.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35116	5.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35117	6.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35118	6.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case
35119	7.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral/Nasal Magill	10/Case

### **HALYARD\* MICROCUFF\* Endotracheal Tube, Oral Curved**



Code	Inner Diameter	Description	Packaging
35161	3.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35162	3.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35163	4.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35164	4.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35165	5.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35166	5.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35167	6.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35168	6.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case
35169	7.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Paediatric Oral Curved	10/Case


NOTE: Please see package inserts for directions for use for each of these products.

# RESPIRATORY HEALTH

## ADULT ENDOTRACHEAL TUBES

### HALYARD\* MICROCUFF\* Oral/Nasal Magill Endotracheal Tube

- Advanced Micro-Thin Polyurethane Cuff Material Designed to Reduce Micro-Aspiration
- Provides a Superior Tracheal Seal and Reduces Leakage Past the Cuff
- Seals at Low Pressure
- Oral/Nasal Magill Curve With Murphy Eye



Code	Inner Diameter	Description	Packaging
35210	5.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35211	5.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35212	6.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35213	6.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35214	7.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35215	7.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35216	8.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35217	8.5mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35218	9.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case
35220	10.0mm	HALYARD* MICROCUFF* Endotracheal Tube, Oral/Nasal Magill, Murphy Eye	10/Case

NOTE: Please see package inserts for directions for use for each of these products.