

RESPIRATORY HEALTH

KIMVENT* BAL CATH*

KIMVENT* BAL CATH* Kit

- Non-Bronchoscopic Bronchoalveolar Lavage Catheter
- Directional tip allows right or left lung sampling
- Maintains PEEP when used with supplied ventilator adapter
- Features a soft, cushioned radiopaque tip for safe sampling from the lower lung



| Code | Outer Diameter | Inner Diameter | Description | Packaging |
|------|----------------|----------------|---|-----------|
| 142 | 16 Fr | 12 Fr | Non-Bronchoscopic Bronchoalveolar Lavage Catheter | 5/Case |
| 143 | 13 Fr | 8 Fr | Non-Bronchoscopic Bronchoalveolar Lavage Catheter | 5/Case |

NOTE: Please see package inserts for directions for use for each of these products.