


INTRODUCER KITS

HALYARD* MIC*, MIC-KEY* Introducer Kits

Kit Includes:

- 1- Gastrointestinal Anchor Set With 4 SAF-T-PEXY* T-Fasteners
 - 1- Haemostat
 - 1- 18 gauge Safety Introducer Needle
 - 1- .035", .9mm J-tip Guidewire
 - 1- #11 Scalpel
 - 1- 12ml Syringe
 - 1- Serial Dilator With Peel-Away Sheath
- Feeding tube not included in kit.




Code	Description	Packaging
98420	MIC* G -12 Introducer Kit Includes 16 Fr Dilator For Use With MIC* 12 Fr Gastrostomy Feeding Tube	1 Kit/Case
98421	MIC* G -14 Introducer Kit Includes 18 Fr Dilator For Use With MIC* 14 Fr Gastrostomy Feeding Tube	1 Kit/Case
98422	MIC* G -16 Introducer Kit Includes 20 Fr Dilator For Use With MIC* 16 Fr Gastrostomy Feeding Tube	1 Kit/Case
98423	MIC* G -18 Introducer Kit Includes 22 Fr Dilator For Use With MIC* 18 Fr Gastrostomy Feeding Tube	1 Kit/Case
98424	MIC* G -20 Introducer Kit Includes 24 Fr Dilator For Use With MIC* 20 Fr Gastrostomy Feeding Tube	1 Kit/Case

HALYARD* MIC* J / TJ Introducer Kits

Kit Includes:

- 1- Gastrointestinal Anchor set with 4 SAF-T-PEXY* T-Fasteners
 - 1- Haemostat
 - 1- 18 gauge Safety Introducer Needle
 - 1- Seeking Catheter 0.38", .97mm Guidewire Compatible
 - 1- #11 Scalpel
 - 1- 12ml Syringe
 - 1- Serial Dilator With Peel-Away Sheath
- Feeding tube not included in kit.



Code	Description	Packaging
98435	MIC* J / TJ - 16 Introducer Kit Includes 20 Fr Dilator For Use With MIC* 16 Fr Jejunal / Transgastric Jejunal Feeding Tube	1 Kit/Case
98436	MIC* J / TJ - 18 Introducer Kit Includes 22 Fr Dilator For Use With MIC* 18 Fr Jejunal / Transgastric Jejunal Feeding Tube	1 Kit/Case

NOTE: Please see package inserts for directions for use for each of these products.

INTRODUCER KITS

HALYARD* MIC* Introducer Kit Components

Code	Description	Packaging
98701	NEW HALYARD* SAF-T-PEXY* T-Fasteners plus Curved Haemostat	3/Pack
98702	HALYARD* 16Fr Serial Dilator, single Recommend using with 12Fr feeding tube	1
98703	HALYARD* 18Fr Serial Dilator, single Recommend using with 14Fr feeding tube	1
98704	HALYARD* 20Fr Serial Dilator, single Recommend using with 16Fr feeding tube	1
98705	HALYARD* 22Fr Serial Dilator, single Recommend using with 18Fr feeding tube	1
98706	HALYARD* 24Fr Serial Dilator, single Recommend using with 20Fr feeding tube	1

NOTE: Please see package inserts for directions for use for each of these products.