



A SMOOTH OPERATION DAY IN, DAY OUT

From your **Trusted Partner** in Acute Care

DAY SURGERY
INFECTION CONTROL
PRODUCT CATALOGUE

AUSTRALIA & NEW ZEALAND

INFECTION PREVENTION SOLUTIONS FOR DAY SURGERY



INFECTION
PREVENTION
AND SURGICAL
SOLUTIONS



EDUCATION &
ADVOCACY



WORKFLOW
AND WASTE
EFFICIENCIES



CLINICAL EXPERTISE
CONSULTANCY















































QUICK LINKS

Infection Prevention in Day Surgery	page 3
STERILISATION CONSUMABLES	page 4
Sterilisation Wrap	page 5
Supplementaries	page 22
Decontamination Gloves	page 24
SURGICAL DRAPES & PACKS	page 25
Abdominal	page 26
Basic Packs	page 28
Basic Drapes	page 29
Breast / Chest	page 31
ENT	page 32
ENT / Dental	page 33
Equipment Covers	page 35
Gynaecology	page 36
Minor Procedure	page 40
Ophthalmic	page 41
Orthopaedic	page 42
Thyroid	page 53
Universal Drapes	page 54
Universal Packs	page 56
Urology	page 59
Supplementaries	page 60
SURGICAL GOWNS	page 62
AAMI Level 4	page 64
AAMI Level 3	page 66
Standard	page 67
CUSTOM PROCEDURE PACKS (CPP)	page 68
PERSONAL PROTECTIVE EQUIPMENT (PPE)	page 70
Gowns	page 72
Masks	page 74
Eyewear	page 79
Gloves	page 80
Headwear	page 85
Footwear	page 86
SURFACE, STAFF & PATIENT CARE	page 87
Surface Care	page 88
Staff & Patient Care	page 92
Patient Care	page 93

INFECTION PREVENTION IN DAY SURGERY

Halyard has a broad range of surgical and infection prevention products to help provide trusted hygienic protection for Patients, Staff, Surfaces and Equipment in Day Surgeries including for:


- Sterilisation / Central Sterile Services Department (CSSD)
- Surgical / Operating Room (O.R.)
- Personal Protective Equipment (PPE)
- Surface Care & Protection
- Fluid Absorption
- Patient Care
- Hand Hygiene

	CSSD / STERILISATION CONSUMABLES	SURGICAL DRAPES & PACKS	SURGICAL GOWNS	PERSONAL PROTECTIVE EQUIPMENT (PPE)	SURFACE CARE	STAFF & PATIENT CARE
HALYARD* STERILIZATION WRAP	 page 5					
HALYARD* CONTAINER FILTERS	 page 22					
HALYARD* TRANSPORT TRAY	 page 22					
HALYARD* ABSORBENT TRAY LINERS	 page 22					
HALYARD TERI* HYDRO ABSORBENT WIPERS	 page 22				 page 89	 page 93
HALYARD* WASHCLOTH EXTRA	 page 23					 page 94
HALYARD TERI* WIPERS	 page 23				 page 89	 page 93
HALYARD ISOWIPE* BACTERICIDAL WIPES	 page 23				 page 90	
HALYARD* DECONTAMINATION GLOVES	 page 24					
HALYARD* SURGICAL DRAPES & PACKS		 page 25				
HALYARD* SURGICAL SUPPLEMENTARIES		 page 60				
HALYARD* O.R. ACCESSORIES		 page 61				
HALYARD* SURGICAL GOWNS			 page 62			
HALYARD* CUSTOM PROCEDURE PACKS	 page 68					
HALYARD* NON-SURGICAL GOWNS	 page 72			 page 72		
HALYARD* SURGICAL MASKS				 page 74		
HALYARD GUARDALL* FACE SHIELD	 page 79			 page 79		
HALYARD SAFEVIEW* GLASSES	 page 79			 page 79		
HALYARD* GLOVES	 page 80			 page 80		
HALYARD* HEADWEAR	 page 85			 page 85		
HALYARD* FOOTWEAR	 page 86			 page 86		
HALYARD* VERSA TOWEL					 page 88	 page 92
HALYARD* UNDERPADS					 page 90	 page 90
HALYARD* BED PROTECTOR					 page 90	 page 90
HALYARD* BEDSHEET					 page 90	
HALYARD* BENCHROLL					 page 91	
HALYARD* SINGLE USE LINEN					 page 91	
HALYARD* NEUTRAL HAND WIPES						 page 92
HALYARD* PROTECTAPAD						 page 93
HALYARD* ICE PACKS						 page 94
HALYARD* I.V. SUPPORT						 page 95

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use. These products are not available for purchase by the general public.



DAY SURGERY PRODUCTS
STERILISATION CONSUMABLES



STERILISATION WRAP

THE HALYARD* STERILIZATION WRAP ADVANTAGE

Since the introduction of single-use sterilisation wrap, health professionals have looked to the market leader for the resources, innovation and high performance solutions you need.

When it comes to your sterilisation packaging, preserving sterility means protecting your patient. HALYARD* Sterilization Wrap with KIMGUARD^ Fabric Technology was the first wrap to receive FDA clearance for all standard sterilisation methods¹, providing you with the reliable solution that you need for optimal patient care and protection.


PROTECTION THROUGH STRENGTH

HALYARD* Sterilization Wrap is made out of polypropylene using a five layer (Spunbond-Meltblown-Meltblown-Meltblown-Spunbond) process; the outer layers of the sterilisation wrap provide exceptional strength and durability.

The inner multi-layers of Meltblown fibres provide an optimal barrier to bacteria and liquids while allowing penetration of gas sterilants (steam, ethylene oxide and plasma). The technical definition of layers differs from that of sheets.

Typical examples of the Spunbond/Meltblown/Spunbond are SMS, SSMMS, and SMMS – abbreviations referring specifically to layer configuration, not number of sheets.

It is not possible to differentiate between the number of meltblown layers; **with everything else equal only the meltblown content (thickness) determines barrier properties – not the number of layers.**



CONFIDENCE


All HALYARD* sterilisation packaging materials are CE-marked and meet AS/NZS 4187 Standards, complying to EN ISO 11607 using tests from EN 868-2.

Next to being strong and being an excellent bacterial barrier, all HALYARD* Sterilization Wrap products have been tested to ensure they meet the Optimum Balance Properties.^{3,4}

YOUR CHECKLIST FOR THE OPTIMUM BALANCE OF PROPERTIES:^{3,4}

Performance Attribute	Clinical Significance
Allows the sterilant to penetrate	• Appropriate sterilisation
Allows air to pass through	
Resists water penetration	• Infection prevention • Barrier efficacy
Resists microbial penetration	
Low linting	
Resists tears, punctures, and abrasion	• Permits aseptic delivery of contents to the sterile field
Does not have excessive "memory"	
Resists ignition / fire	• Fire safety

LOOK
FOR THE
DAISY
PATTERN



STERILISATION WRAP

PROTECTION THROUGH AN OPTIMAL BACTERIAL BARRIER

HALYARD* sterilisation packaging material is treated with unique POWERGUARD* Technology. The Halyard process infuses the wrap fabric with unique electric micro-fields surrounding the Meltblown fibres.

This varied network of charged fibres forms gradients within the fabric structure, enabling it to be highly effective at microbial capture and retention – so essential in preventing fibre release of captured contaminants over time.

The POWERGUARD* treatment is fabric neutral (pluses and minuses staying separated) and maintains these high capture and retention properties through sterilisation, manipulation and storage.



Reusable sterilisation wrap



SMS sterilisation wrap




HALYARD* ONE-STEP* Sterilization Wrap H400²

These images illustrate the effect of POWERGUARD* on sterile barrier materials after BFE-testing (ASTM F2101): the plates have been covered respectively with a reusable sterilisation wrap, two layers of SMS sterilisation wrap, and ONE-STEP* Sterilization Wrap H400. The fabric covered plates were then challenged with *Staphylococcus aureus* bacteria. After twenty-four hours, the wrap is removed and the CFUs are counted to determine the number of bacteria that penetrated through the different sterilisation wrap fabrics. Results are reported in the number of bacteria that are prevented from passing through the material.

BARRIER INTEGRITY

Percent Filtration Efficiency (BFE)	Number Bacteria Passed Through
50%	1,100
70%	660
90%	220
95%	110
96%	88
97%	66
98%	44
99%	22
≥ 99.9%[†]	0–2

[†] 100% BFE reporting is not permitted per good statistical practices



STERILISATION WRAP

CLINICALLY SUPERIOR PRODUCTS

1st and only wrap manufacturer with FDA clearance for ALL approved sterilisation methods.

Only Halyard has 510(k) clearance from the FDA for all standard modalities including Pre-vacuum Steam, STERRAD®, STERIS®, V-PRO, EO and Gravity. We also include simplified *Instructions for Use* in every case of wrap.

We continually re-set the bar for product quality and performance through clinical and scientific research directed by on-staff medical professionals as well as peer-to-peer consultation.

ONE STEP* QUICK CHECK* CLEARANCES⁵

PRE-VACUUM STEAM
ETHYLENE OXIDE (EO)
GRAVITY STEAM
STERIS® V-PRO® 1
STERIS® V-PRO® 1 PLUS
STERIS® V-PRO® MAX
STERIS® V-PRO® 60
STERRAD® 100NX
STERRAD® NX
STERRAD® 100S
STERRAD® 200
STERRAD® 50
STERILUCENT

SMART-FOLD* CLEARANCES⁶

PRE-VACUUM STEAM
ETHYLENE OXIDE (EO)
STERRAD® 100NX
STERRAD® NX
STERRAD® 100S
STERIS® V-PRO® 60
STERIS® V-PRO® 1
STERIS® V-PRO® 1 PLUS
STERIS® V-PRO® MAX
STERIS® V-PRO® MAX 2

DEEP EXPERTISE AND SERVICE



Our Excellence in Sterile Processing (ESP) curriculum offers more than 25 courses that address Central Sterilisation Processing Department requirements.




To help drive safety and efficiency, we also offer Sterilisation Packaging Utilisation Reviews (SPURS).

Our team is available to assist you, at no additional cost to your facility, with:

- **In-service training**
- **Evaluations**
- **Clinical Consultancy**
 - Tear Identification and Handling Technique
 - Wet Packs
 - Event Related Sterility Maintenance
- **Technical Consultation Services**
 - Wet Packs Audit
 - Wrap Tear Audit

- **Commitment to Excellence and Sustainability**

As a company that both manufactures and sells our own product, we strive to deliver the highest quality and sustainability standards for our customers. We do this through sustainable sourcing, working with our customers to reduce waste through effective product use and purchasing, and supporting responsible disposal and recycling.



LOOK
FOR THE
DAISY
PATTERN

STERILISATION WRAP


WHY DOUBLE-WRAPPING IS IMPERATIVE FOR OPTIMAL BARRIER PROTECTION

"Sterility assurance levels of 10 to the minus sixth (one part per million) is the level of microorganism colony growth that would be considered not harmful to humans or able to produce disease" ¹³

Only double-wrapping provides these assurance levels.¹⁴


Due to the porosity requirements of SMS for sterilant penetration, no barrier material can be 100% effective in filtering bacteria, therefore assuming 99.9% filtration efficiency there is a:

One in a thousand chance that items wrapped in a single sheet of wrap can become contaminated.




Based on 100 trays processed each day, the potential for contamination is 1 event every 10 days.

One in a million chance that items wrapped in two sheets of wrap can become contaminated.




Based on 100 trays processed each day, potential for contamination is 1 event every 27 years.

LOOK
FOR THE
DAISY
PATTERN




STERILISATION WRAP

WHY BONDED WRAP EXHIBITS THE SAME IF NOT BETTER BENEFITS THAN SEQUENTIAL WRAPPING



ONE-STEP* and QUICK CHECK* Sterilization Wrap allow for single-step wrapping. Studies conducted by Kurt-Salmon Associates demonstrate 49% time saving when wrapping and 48% time saving upon opening.⁷




Barrier protection is equal to packaging with 2 single sheets as proven by TNO.^{15,16}

Objective: To objectively evaluate the barrier properties of an assembled, sealed and sterilised pack TNO (The Netherlands Organization for Applied Scientific Research) studied the barrier properties of a selection of sterilisation wrapping materials.

Results: Results are expressed as % retention with zero being no barrier and a retention greater than 99.9% being an efficient barrier. All weights of HALYARD* Sterilization Wrap (both ONE-STEP* and two sheets of HALYARD* Sequential Sterilization Wrap) exhibit the same excellent barrier characteristics.

"...using two barrier-type wrappers provided a tortuous pathway to impede microbial migration and permits ease of presentation to the sterile field without compromising sterility."¹⁷



STERILISATION WRAP

At Halyard, we believe that there should be a sterilisation wrap available for every need and for every user. Within the HALYARD* sterilisation portfolio we offer:



HALYARD* SEQUENTIAL STERILIZATION WRAP

PROTECT EVERY TRAY WITH CONFIDENCE – SEQUENTIAL WRAPPING

Extensive independent laboratory testing assures that all HALYARD* Sterilization Wrap delivers the **optimum balance of properties**³ needed to ensure the sterility of your instruments:

- Strength and resistance to tearing³
- Resistance to linting and abrasion³
- Effective microbial barrier³
- Effective sterilant penetration³
- Post sterilisation testing for confidence in sterility of instruments with Event Related Sterility Maintenance Studies.³

Available in four grades: H200, H300, H500 and H600; the higher the number, the stronger the wrap.



QUICK CHECK* STERILIZATION WRAP

WRAP WITH CONFIDENCE YOU CAN SEE – SIMULTANEOUS WRAPPING

QUICK CHECK* Sterilization Wrap has the added benefit of making inspecting trays easier and faster⁷.

- Delivers the **optimum balance of properties**³
- Two sheets of wrap thermally bonded on both sides allowing you to save up to 49% of time⁷
- Each sheet is a contrasting colour – one white and one blue⁸
- Post sterilisation testing for confidence in sterility of instruments with Final Pack⁹ and Maintenance of Package Integrity¹⁰ tests.

Available in three grades: H200, H400 and H600; the higher the number, the stronger the wrap.



SMART-FOLD* STERILIZATION WRAP


CONFIDENCE IN THE MOST CHALLENGING TRAYS – STRONGER, FASTER, EASIER AND MORE CONSISTENT WRAP BY DESIGN

SMART-FOLD* Sterilization Wrap delivers benefits beyond QUICK CHECK* Sterilization Wrap:


- Delivers the **optimum balance of properties**⁵
- Patented design with triple strength reinforcement zones to protect your heaviest trays and loaner sets
- 64% more durable¹¹
- 86% faster inspection time¹²
- 59% faster wrapping time¹²
- Post sterilisation testing for confidence in sterility of instruments with Final Pack⁹ and Maintenance of Package Integrity¹⁰ tests.

Available in two grades: H450 and H650; the higher the number, the stronger the wrap.


All HALYARD* sterilisation products exist in different sizes and strengths depending on the weight and size of the materials to be packaged:




H200




H300



H400





H500



H600

LOOK
FOR THE
DAISY
PATTERN





HALYARD* SEQUENTIAL STERILIZATION WRAP



- Single sheet polypropylene:
- Features POWERGUARD* Technology for highest bacterial filtration efficiency
- Low in lint
- Different sizes and grades depending on the weight and size of the materials to be packaged
- Water resistant
- Tested for event-related and time-related sterility maintenance.
- Flammability resistant

- Validated with the most current sterilisation modalities:
 - Pre-Vacuum Steam
 - Gravity Steam Sterilisation
 - Ethylene Oxide (EO)
 - STERRAD® Sterilisation Systems†
 - STERILUCENT
 - STERIS® Sterilisation Systems††

† STERRAD® NX, STERRAD® 100NX, STERRAD® 100S, STERRAD® 200, STERRAD® 50.

†† STERIS® V-PRO® 1, STERIS® V-PRO® 1 PLUS, STERIS® V-PRO® MAX, STERIS® V-PRO® 60.

HALYARD* SEQUENTIAL STERILIZATION WRAP H200

Code	Size	Packaging
68012	30cm x 30cm	Units/Pack = 500 Packs/Case = 2 Units/Case = 1000
68018	45cm x 45cm	Units/Pack = 500 Packs/Case = 2 Units/Case = 1000

HALYARD* SEQUENTIAL STERILIZATION WRAP H300

Code	Size	Packaging
10824	60cm x 60cm	Units/Pack = 250 Packs/Case = 2 Units/Case = 500
10830	76cm x 76cm	Units/Pack = 100 Packs/Case = 3 Units/Case = 300
10836	91cm x 91cm	Units/Pack = 100 Packs/Case = 3 Units/Case = 300
10840	101cm x 101cm	Units/Pack = 50 Packs/Case = 5 Units/Case = 250
10848	121cm x 121cm	Units/Pack = 50 Packs/Case = 2 Units/Case = 100
10854	137cm x 137cm	Units/Pack = 50 Packs/Case = 2 Units/Case = 100
10872	137cm x 182cm	Units/Pack = 50 Packs/Case = 2 Units/Case = 100



HALYARD* SEQUENTIAL STERILIZATION WRAP H500

Code	Size	Packaging
68130	76cm x 76cm	Units/Pack = 125 Packs/Case = 2 Units/Case = 250
68136	91cm x 91cm	Units/Pack = 75 Packs/Case = 2 Units/Case = 150
68148	121cm x 121cm	Units/Pack = 50 Packs/Case = 1 Units/Case = 50
68154	137cm x 137cm	Units/Pack = 50 Packs/Case = 1 Units/Case = 50
68172	137cm x 182cm	Units/Pack = 50 Packs/Case = 1 Units/Case = 50
68160	152cm x 152cm	Units/Pack = 25 Packs/Case = 2 Units/Case = 50



HALYARD* SEQUENTIAL STERILIZATION WRAP H600

Code	Size	Packaging
68648	121cm x 121cm	Units/Pack = 48 Packs/Case = 1 Units/Case = 48
68672	137cm x 182cm	Units/Pack = 24 Packs/Case = 2 Units/Case = 48



QUICK CHECK* STERILIZATION WRAP

- Two sheets of polypropylene wrap bonded along 2 sides
- Easy breach identification – each sheet contrasting colour, one white and one blue
- Features POWERGUARD* technology for highest bacterial filtration efficiency
- All the double sheet ease and time-savings of ONE-STEP* Sterilization Wrap
- Low in lint
- Different sizes and grades depending on the weight and size of the materials to be packaged
- Water resistant
- Tested for event-related and time-related sterility maintenance
- Flammability resistant



- Validated with the most current sterilisation modalities:
 - Pre-Vacuum Steam
 - Gravity Steam Sterilisation
 - Ethylene Oxide (EO)
 - STERRAD® Sterilisation Systems†
 - STERILUCENT
 - STERIS® Sterilisation Systems††

† STERRAD® NX, STERRAD® 100NX, STERRAD® 100S, STERRAD® 200, STERRAD® 50.

†† STERIS® V-PRO® 1, STERIS® V-PRO® 1 PLUS, STERIS® V-PRO® MAX, STERIS® V-PRO® 60.

QUICK CHECK* STERILIZATION WRAP H200

Code	Size	Packaging
34167	30cm x 30cm	Applications/Pack = 240 Packs/Case = 2 Applications/Case = 480
34193	45cm x 45cm	Applications/Pack = 120 Packs/Case = 2 Applications/Case = 240
34174	60cm x 60cm	Applications/Pack = 120 Packs/Case = 2 Applications/Case = 240
34175	76cm x 76cm	Applications/Pack = 48 Packs/Case = 3 Applications/Case = 144
34185	91cm x 91cm	Applications/Pack = 48 Packs/Case = 3 Applications/Case = 144
34186	101cm x 101cm	Applications/Pack = 24 Packs/Case = 5 Applications/Case = 120
34149	121cm x 121cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34187	137cm x 137cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34168	137cm x 182cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48

QUICK CHECK* Sterilization Wrap continues over the page...

QUICK CHECK* STERILIZATION WRAP H400



Code	Size	Packaging
34183	60cm x 60cm	Applications/Pack = 60 Packs/Case = 2 Applications/Case = 120
34158	91cm x 91cm	Applications/Pack = 36 Packs/Case = 2 Applications/Case = 72
34184	101cm x 101cm	Applications/Pack = 24 Packs/Case = 3 Applications/Case = 72
34157	114cm x 114cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34145 [†]	121cm x 121cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34146 [†]	137cm x 137cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34191	137cm x 182cm	Applications/Pack = 12 Packs/Case = 2 Applications/Case = 24

QUICK CHECK* STERILIZATION WRAP H500



Code	Size	Packaging
34165	76cm x 76cm	Applications/Pack = 60 Packs/Case = 2 Applications/Case = 120

QUICK CHECK* STERILIZATION WRAP H600



Code	Size	Packaging
34151	91cm x 91cm	Applications/Pack = 36 Packs/Case = 2 Applications/Case = 72
34163	101cm x 101cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34164	114cm x 114cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34147 [†]	121cm x 121cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34142	137cm x 137cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34197	137cm x 182cm	Applications/Pack = 12 Packs/Case = 2 Applications/Case = 24

NOTE: Instructions for use are available at www.halyardhealth.com.au
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.



SMART-FOLD* STERILIZATION WRAP

- Unique shape and revolutionary design provides consistent wrapping every time
- White inner sheet aids in breach identification
- Features POWERGUARD* technology for highest bacterial filtration efficiency
- Triple thickness reinforcement zones to protect edges and corners
- Low in lint
- Different sizes and grades depending on the weight and size of the materials to be packaged
- Water resistant
- Tested for event-related and time-related sterility maintenance
- Flammability resistant
- 59% faster wrapping compared to ONE-STEP* Wrap†
- 86% faster inspections†
- 22% less bulk††

- Validated with the most current sterilisation modalities:
 - Pre-Vacuum Steam
 - Ethylene Oxide (EO)
 - STERRAD® Sterilisation Systems†
 - STERILUCENT

† STERRAD® NX, STERRAD® 100NX, STERRAD® 100S, STERRAD® V-PRO® 60, STERRAD® V-PRO® 1, STERRAD® V-PRO® 1 PLUS, STERRAD® V-PRO® MAX, STERRAD® V-PRO® MAX 2.



SMART-FOLD* STERILIZATION WRAP H450

Code	Size (includes side-to-side and pull tabs)	Packaging
14268	Size 1: 101cm x 119cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48
14275	Size 3: 101cm x 139cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48
14279	Size 4: 71cm x 116cm	Units/Pack = 12 Packs/Case = 6 Units/Case = 72
14324	Size 8: 55cm x 114cm	Units/Pack = 24 Packs/Case = 3 Units/Case = 72
14313	Size 10: 121cm x 154cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48

SMART-FOLD* STERILIZATION WRAP H650

Code	Size (includes side-to-side and pull tabs)	Packaging
14271	Size 1: 101cm x 119cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48
14277	Size 3: 101cm x 139cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48
14281	Size 4: 71cm x 116cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48

† Comparison of Smart-Fold* Wrap to One-Step* Wrap and Sequential Wrap: Kurt Salmon & Associates, Time and Motion Studies on the Labor Requirements of the Smart-Fold, One-Step, and Sequential Sterilization Wraps, June 2011.

†† Halyard product specifications.

NOTE: Instructions for use are available at www.halyardhealth.com.au

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

At Halyard we understand that your most challenging trays are usually your heaviest, most complex orthopaedic trays. As these trays have the largest volume, you need a range of wraps that have been designed to withstand the weight and deliver confidence in the OR. These wraps are not just limited to orthopaedic – this range would suit any heavy tray.



QUICK CHECK* STERILIZATION WRAP H600

Code	Size	Packaging
34151	91cm x 91cm	Applications/Pack = 36 Packs/Case = 2 Applications/Case = 72
34163	101cm x 101cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34164	114cm x 114cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34147 [†]	121cm x 121cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34142	137cm x 137cm	Applications/Pack = 24 Packs/Case = 2 Applications/Case = 48
34197	137cm x 182cm	Applications/Pack = 12 Packs/Case = 2 Applications/Case = 24

SMART-FOLD* STERILIZATION WRAP H450

Code	Size (includes side-to-side and pull tabs)	Packaging
14268	Size 1: 101cm x 119cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48
14275	Size 3: 101cm x 139cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48
14279	Size 4: 71cm x 116cm	Units/Pack = 12 Packs/Case = 6 Units/Case = 72
14324	Size 8: 55cm x 114cm	Units/Pack = 24 Packs/Case = 3 Units/Case = 72
14313	Size 10: 121cm x 154cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48

SMART-FOLD* STERILIZATION WRAP H650

Code	Size (includes side-to-side and pull tabs)	Packaging
14271	Size 1: 101cm x 119cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48
14277	Size 3: 101cm x 139cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48
14281	Size 4: 71cm x 116cm	Units/Pack = 12 Packs/Case = 4 Units/Case = 48

STERILISATION POUCHES

From the brand you know and trust in sterilisation wrap, Halyard now offers that same confidence in performance through our extensive range of reliable, high quality HALYARD* Sterilization Pouches to meet your sterility needs.

- Puncture and abrasion resistance with the specially formulated HALYARD* multilayer transparent film¹⁸
- Quick check of seal integrity with colour changing film and contrasting indicators¹⁸
- Triple rail heat seals¹⁸
- Exceptional clean and fibre-free peel¹⁸
- Easy aseptic presentation and opening with a top thumb notch
- Heat-sealed corners on pre-formed pouches
- Tested for event-related and time-related sterility maintenance¹⁸
- All chemical indicators and printing external to the packing area¹⁸
- For Steam and Ethylene Oxide (EO) sterilisation, use Paper/Film pouches¹⁸:
 - 68gsm premium medical grade paper
 - Indicator changes from green to purple for STEAM sterilisation:
- Indicator changes from yellow to brown for EO sterilisation:
- For Low Temperature sterilisation with Hydrogen Peroxide (H₂O₂), use Tyvek®/Film pouches:¹⁸⁻²⁰
 - Tyvek® 1073B material
 - Indicator changes from blue to pink:

GREEN

>

PURPLE

YELLOW

>

BROWN

BLUE

>

PINK





HALYARD* HEAT-SEAL PAPER/FILM STERILIZATION POUCHES

Code	Size	Compatible Sterilisation Modality	Packaging
48516	50mm x 250mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48517	75mm x 150mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48518	75mm x 200mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48519	75mm x 250mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48496	100mm x 200mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48497	100mm x 250mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48498	100mm x 270mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48499	100mm x 300mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48506	150mm x 300mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48507	150mm x 380mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48508	150mm x 400mm	STEAM / EO	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48510	205mm x 400mm	STEAM / EO	Units/Box = 250 Boxes/Case = 4 Units/Case = 1000
48512	250mm x 380mm	STEAM / EO	Units/Box = 250 Boxes/Case = 4 Units/Case = 1000
48513	250mm x 500mm	STEAM / EO	Units/Box = 250 Boxes/Case = 2 Units/Case = 500

HALYARD* HEAT-SEAL PAPER/FILM STERILIZATION REELS

Code	Size	Compatible Sterilisation Modality	Packaging
48525	50mm x 200m	STEAM / EO	Reels/Case = 2
48526	75mm x 200m	STEAM / EO	Reels/Case = 2
48521	100mm x 200m	STEAM / EO	Reels/Case = 1
48522	150mm x 200m	STEAM / EO	Reels/Case = 1
48523	200mm x 200m	STEAM / EO	Reels/Case = 1
48524	250mm x 200m	STEAM / EO	Reels/Case = 1


HALYARD* HEAT-SEAL PAPER/FILM GUSSETED STERILIZATION POUCHES

Code	Size	Compatible Sterilisation Modality	Packaging
48493	200mm x 55mm x 400mm	STEAM / EO	Units/Box = 100 Boxes/Case = 10 Units/Case = 1000
48494	150mm x 50mm x 100m	STEAM / EO	Units/Case = 2



HALYARD* HEAT-SEAL TYVEK®/FILM STERILIZATION POUCHES

Code	Size	Compatible Sterilisation Modality	Packaging
48532	75mm x 210mm	PLASMA (H ₂ O ₂)	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48527	100mm x 260mm	PLASMA (H ₂ O ₂)	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48528	150mm x 320mm	PLASMA (H ₂ O ₂)	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48529	200mm x 400mm	PLASMA (H ₂ O ₂)	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48530	200mm x 500mm	PLASMA (H ₂ O ₂)	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000
48531	200mm x 650mm	PLASMA (H ₂ O ₂)	Units/Box = 250 Boxes/Case = 8 Units/Case = 2000



HALYARD* HEAT-SEAL TYVEK®/FILM STERILIZATION REELS

Code	Size	Compatible Sterilisation Modality	Packaging
48534	150mm x 100m	PLASMA (H ₂ O ₂)	Reels/Box = 1 Boxes/Case = 8 Reels/Case = 8
48535	200mm x 100m	PLASMA (H ₂ O ₂)	Reels/Box = 1 Boxes/Case = 4 Reels/Case = 4
48536	250mm x 100m	PLASMA (H ₂ O ₂)	Reels/Box = 1 Boxes/Case = 4 Reels/Case = 4

HALYARD* SELF-SEAL PAPER/FILM STERILIZATION POUCHES


Code	Size	Compatible Sterilisation Modality	Packaging
48548	90mm x 130mm	STEAM / EO	Units/Box = 200 Boxes/Case = 10 Units/Case = 2000
48549	90mm x 200mm	STEAM / EO	Units/Box = 200 Boxes/Case = 10 Units/Case = 2000
48550	90mm x 250mm	STEAM / EO	Units/Box = 200 Boxes/Case = 10 Units/Case = 2000
48541	200mm x 350mm	STEAM / EO	Units/Box = 200 Boxes/Case = 5 Units/Case = 1000
48544	250mm x 400mm	STEAM / EO	Units/Box = 200 Boxes/Case = 5 Units/Case = 1000
48545	300mm x 450mm	STEAM / EO	Units/Box = 200 Boxes/Case = 5 Units/Case = 1000

HALYARD* SELF-SEAL TYVEK®/FILM STERILIZATION POUCHES

Code	Size	Compatible Sterilisation Modality	Packaging
48552	150mm x 320mm	PLASMA (H ₂ O ₂)	Units/Box = 200 Boxes/Case = 10 Units/Case = 2000
48554	75mm x 200mm	PLASMA (H ₂ O ₂)	Units/Box = 200 Boxes/Case = 10 Units/Case = 2000

STERILISATION CONSUMABLES

SUPPLEMENTARIES



HALYARD* STERILIZATION WRAP CONTAINER FILTERS

- H200 Sterilisation Wrap
- 100% Polypropylene with SMS Structure
- Latex-Free
- Anti-Static Treated

Code	Description	Packaging
68507	HALYARD* Container Filter, Aesculap*, 19cm Round	Eaches/Case = 1,000
68509	HALYARD* Container Filter, Genesis*, 22.9cm Square	Eaches/Case = 1,000



HALYARD* TRANSPORT TRAY

- Polypropylene with SMS Structure
- Blue

Code	Description	Packaging
KL-1500	Transport Tray, Universal, 29cm x 55cm	Eaches/Case = 10
KL-1500XL	Transport Tray, Extra Large, 67.6cm x 41.9cm	Eaches/Case = 5



HALYARD* ABSORBENT TRAY LINERS

- High Absorbency – absorbs up to six times its own weight in water²³ to help wick away excessive moisture during steam sterilisation
- Low linting²⁶
- Cut to fit most standard Sterilisation Trays – may be folded to accommodate different tray sizes
- ARTG 332698

Code	Size	Description	Packaging
4330	21.5cm x 24cm	HALYARD* Absorbent Tray Liners Small	Eaches/Pack = 124 Packs/Case = 4 Eaches/Case = 496
4331	24cm x 36.5cm	HALYARD* Absorbent Tray Liners Medium	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4332	30cm x 50.5cm	HALYARD* Absorbent Tray Liners Large	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450



HALYARD TERI* HYDRO ABSORBENT WIPERS

- High Absorbency – absorbs up to six times its own weight in water²³ to help wick away excessive moisture during steam sterilisation
- Ultra low linting²⁶
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698

Code	Size	Description	Packaging
4476	32.5cm x 33cm	TERI* Hydro Absorbent Wipers Small	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4478	32.5cm x 60cm	TERI* Hydro Absorbent Wipers Large	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450

NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.

STERILISATION CONSUMABLES

SUPPLEMENTARIES



HALYARD* WASHCLOTH EXTRA

- Soft, thick and super absorbent - absorbs up to eight times its weight in water²³
- Strong and durable enough to not tear easily²² in wet or dry applications
- Low linting²⁶
- Clear resealable bag to help prevent dust and other particles from contaminating cloths

Code	Size	Description	Packaging
2711	35cm x 40cm	HALYARD* Washcloth Extra	Eaches/Pack = 20 Packs/Case = 30 Eaches/Case = 600



TERI* WIPERS

- Soft 4 ply wiper with tough polyester scrim that provide excellent absorbency and strength even when wet²²
- Thick 4 ply tissue provides a substantial barrier between your hands and contaminants
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698

Code	Size	Description	Packaging
4465	31.5cm x 34cm	TERI* Wipers Small	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600
4467	31.5cm x 34cm	TERI* Wipers Small Bulk Pack	Eaches/Pack = 450
4466	31.5cm x 61cm	TERI* Wipers Large	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600

HALYARD ISOWIPE* BACTERICIDAL WIPES

- Contains 70% Isopropyl Alcohol to provide hospital grade disinfection of clean surfaces
- ISOWIPE* Bactericidal Wipe is entered in the ARTG, 342811
- Kills COVID-19 and SARS virus in 60 seconds
- Kills Bacteria (S. aureus, E. coli, Ps. Aeruginosa, Pr. vulgaris)
- Kills Fungi (Yeast and Mould) in 10 minutes

Always read the label and follow the directions for use.

Code	Size	Description	Packaging
6835	42cm x 14cm	ISOWIPE* Bactericidal Wipes Canister	Eaches/Canister = 75 Canisters/Case = 12 Eaches/Case = 900
6836	42cm x 14cm	ISOWIPE* Refill Pack	Eaches/Refill = 75 Refills/Case = 12 Eaches/Case = 900
6837	21cm x 14cm	ISOWIPE* Minis Canister	Eaches/Canister = 150 Canisters/Case = 12 Eaches/Case = 1800
6838	21cm x 14cm	ISOWIPE* Minis Refill Pack	Eaches/Refill = 150 Refills/Case = 12 Eaches/Case = 1800



4989	17.8cm H x 9.4cm D	Stainless Steel Wall Bracket	1
------	--------------------	------------------------------	---

NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.

STERILISATION CONSUMABLES

DECONTAMINATION GLOVES



HALYARD* PURPLE NITRILE* MAX EXAM GLOVES

- Ideal for CSSD
- Purple, Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Longer length: 16" (406mm)
- CE Marked
- Enhanced Wet and Dry Grip
- Heavier weight: 20 grams
- Increased fingertip thickness
- Tested Against 52 Chemotherapy and other Hazardous Agents




Code	Description	Packaging
44992	PURPLE NITRILE* MAX Exam Gloves, Small	Eaches/Box = 50 Boxes/Case = 8
44993	PURPLE NITRILE* MAX Exam Gloves, Medium	Eaches/Box = 50 Boxes/Case = 8
44994	PURPLE NITRILE* MAX Exam Gloves, Large	Eaches/Box = 50 Boxes/Case = 8
44995	PURPLE NITRILE* MAX Exam Gloves, X-Large	Eaches/Box = 50 Boxes/Case = 8




DAY SURGERY PRODUCTS

SURGICAL DRAPES & PACKS




LAPAROSCOPY / LAPAROTOMY DRAPES



LAPAROSCOPIC / CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 30cm x 34cm Adjustable Fenestration, Lateral Offset
- Hook and Loop Tube and Line Control
- CONTROL* Instrument Pad Reinforcement
- Armboard Covers


Code	Description	Packaging
89229	Laparoscopic / Cholecystectomy Drape, Sterile	Items/Case = 12



LAPAROSCOPIC ABDOMINAL DRAPE, 254cm x 193cm x 315cm

- 30cm x 33cm Abdominal Fenestration
- CONTROL* Fabric Reinforcement
- Clear Instrument and Line Control Pockets
- Hook and Loop Tube and Line Control
- Armboard Covers

Code	Description	Packaging
89233	Laparoscopic Abdominal Drape, Sterile	Items/Case = 15
79233	Laparoscopic Abdominal Drape, Non-Sterile	Items/Case = 18



LAPAROSCOPIC DRAPE, 264cm x 89cm x 315cm

- 31cm x 33cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Clear Film Pouch
- Tube Holders
- Armboard Covers
- Under Leg Wings

Code	Description	Packaging
29240	Laparoscopic Drape, Sterile	Items/Case = 15
29241	Laparoscopic Drape, Non-Sterile	Items/Case = 19

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


LAPAROSCOPY / LAPAROTOMY DRAPES



HALYARD BASICS* LAPAROSCOPIC/CHOLECYSTECTOMY DRAPE, 264cm x 193cm x 305cm

- 31cm x 33cm Abdominal Fenestration
- Hook and Loop Tube and Line Control
- Fabric Reinforcement
- Armboard Covers


Code	Description	Packaging
89538	Laparoscopic/Cholecystectomy Drape, Sterile	Items/Case = 25
79538	Laparoscopic/Cholecystectomy Drape, Non-Sterile	Items/Case = 26



PAEDIATRIC LAPAROTOMY DRAPE, 152cm x 274cm

- 7cm x 20cm Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89241	Paediatric Laparotomy Drape, Sterile	Items/Case = 15
79240	Paediatric Laparotomy Drape, Non-Sterile	Items/Case = 25



LAPAROTOMY DRAPE WITH INCISE, 254cm x 183cm x 305cm

- 20cm x 33cm Adhesive Incise Area
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89282	Laparotomy Drape with Incise, Sterile	Items/Case = 15
79282	Laparotomy Drape with Incise, Non-Sterile	Items/Case = 15

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

BASIC PACK VI

Code	Description	Packaging	
88161	Basic Pack VI	Items/Case = 7	
Sterile Component Code	Non-Sterile Component Code	Contents	
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
95111	71311	1-ULTRA Surgical Gown, In Overwrap with Towel, Large	
95121	71321	1-ULTRA Surgical Gown, with Towel, X-Large	
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm	
89731	79731	4-Utility Drapes, 66cm x 38cm, with Tape	
		1-Suture Bag	
47614	79101	1-Half Drape, 147cm x 101cm	
47622	79121	1-Large Drape, 139cm x 193cm	

BASIC PACK XV

Code	Description	Packaging	
88178	Basic Pack XV	Items/Case = 15	
Sterile Component Code	Non-Sterile Component Code	Contents	
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm	
		1-Suture Bag	
47622	79121	1-Large Drape, 139cm x 193cm	

SET-UP PACK I

Code	Description	Packaging	
88171	Set-Up Pack I	Items/Case = 14	
Sterile Component Code	Non-Sterile Component Code	Contents	
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm	
95111	71311	1-Ultra Surgical Gown In Overwrap with Towel, Large	
95121	71321	1-ULTRA Surgical Gown, X-Large	
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm	
		1-Suture Bag	

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.




UTILITY DRAPE, 90cm x 90cm

Code	Description	Packaging
89353	Utility Drape, Sterile	Items/Case = 120
76323	Utility Drape, Non-Sterile	Items/Case = 200




HALF DRAPE, 147cm x 101cm

Code	Description	Packaging
47614	Half Drape, Sterile	Items/Case = 60
47615	Half Drape, Non-Sterile	Items/Case = 1000




MEDIUM DRAPE, 193cm x 112cm

Code	Description	Packaging
89111	Medium Drape, Sterile	Items/Case = 5




MEDIUM DRAPE WITH TAPE, 193cm x 105cm

Code	Description	Packaging
89023	Medium Drape with Tape, Sterile	Items/Case = 50



MEDIUM DRAPE, 185cm x 101cm

Code	Description	Packaging
47619	Medium Drape, Non-Sterile	Items/Case = 800



LARGE DRAPE, 193cm x 139cm

Code	Description	Packaging
47622	Large Drape, Sterile	Items/Case = 32
47624	Large Drape, Non-Sterile	Items/Case = 600

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.



THREE-QUARTER DRAPE, 152cm x 193cm


- CONTROL* Plus Fabric Reinforcement
- No Adhesive Tape

Code	Description	Packaging
79147	Three-Quarter Drape, Non-Sterile	Items/Case = 300




THREE-QUARTER DRAPE, 139cm x 193cm

Code	Description	Packaging
47628	Three-Quarter Drape, Sterile	Items/Case = 22



X-LARGE DRAPE, 193cm x 243cm


Code	Description	Packaging
47630	X-Large Drape, Sterile	Items/Case = 24
47631	X-Large Drape, Non-Sterile	Items/Case = 32



SMALL ABSORBENT DRAPE WITH TAPE, 57cm x 60cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89828	Small Absorbent Drape with Tape, Sterile	Items/Case = 72
79828	Small Absorbent Drape with Tape, Non-Sterile	Items/Case = 108



STANDARD ABSORBENT DRAPE WITH TAPE, 91cm x 91cm

- CONTROL* Plus Absorbent Fabric

Code	Description	Packaging
89343	Standard Absorbent Drape with Tape, Sterile	Items/Case = 50
79343	Standard Absorbent Drape with Tape, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.



**BREAST DRAPE,
254cm x 193cm x 315cm**

- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Armboard Covers

Code	Description	Packaging
89236	Breast Drape, Sterile	Items/Case = 10




**CHEST U-DRAPE,
193cm x 305cm**

- CONTROL* Plus Fabric Reinforcement
- 30cm x 76cm Split
- Tube Holders


Code	Description	Packaging
79335	Chest U-Drape, Non-Sterile	Items/Case = 22

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.




EENT PACK I

Code	Description	Packaging
88711	EENT Pack I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, with Towel, Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
		1-Absorbent Towel, 38cm x 56cm
		1-Suture Bag
89431	79430	1-Head Drape, 127cm x 91cm
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split



EENT PACK II

Code	Description	Packaging
88716	EENT Pack II	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown In Overwrap, with Towel, Large
95111	71311	1-ULTRA Surgical Gown with Towel, Large
95211	74911	1-ULTRA Surgical Gown, Fabric Reinforced, with Towel, Large
		1-Microsurgery Mayo Stand Cover
		1-Suture Bag
89431	79430	1-Head Drape, 127cm x 91cm
89352	79351	1-EENT Split Drape with CONTROL* Reinforcement, 193cm x 315cm with 7cm x 71cm Split



ENT PACK III

Code	Description	Packaging
88223	ENT Pack III	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
89603	79603	1-Back Table Cover, Zone-Reinforced, 152cm x 191cm
		1-Mayo Stand Cover 80cm x 145cm
		1-Suture Bag
89431	79430	1-Head Drape 91cm x 127cm:
		• Turban
		• Positioning Tape
89352	79351	1-ENT Split Drape 193cm x 315cm:
		• with 7 x 71cm Split
		• Positioning Tape
		• CONTROL* Fabric Reinforcement
		• Tube Holders

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.



EENT PACK IV


Code	Description	Packaging
88719	EENT Pack IV	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Surgical Gown with Set-in sleeves, In Overwrap, with Towel, Large
89701	79720	1-Surgical Gown, with Reinforced Sleeves, Large
89601	79601	1-Absorbent Towel, 38cm x 56cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89438	79438	1-Head Drape, 104cm x 81cm 1-EENT Split Drape with Instrument Pad Reinforcement, 193cm x 315cm with 8cm x 64cm Split



EENT SPLIT DRAPE, 193cm x 315cm

- CONTROL* Fabric Reinforcement
- 7cm x 71cm Split
- Tube Holders
- Instrument Pad


Code	Description	Packaging
89352	EENT Split Drape, Sterile	Items/Case = 18
79351	EENT Split Drape, Non-Sterile	Items/Case = 25



HEAD DRAPE, 127cm x 91cm

- Turban, 83cm x 50cm

Code	Description	Packaging
89431	Head Drape, Sterile	Items/Case = 38
79430	Head Drape, Non-Sterile	Items/Case = 60



HEAD DRAPE, 104cm x 81cm

- Detachable Turban, 91cm x 61cm

Code	Description	Packaging
89438	Head Drape, Sterile	Items/Case = 38

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.



DENTAL DRAPE, 120cm x 150cm

- Teardrop Fenestration 10cm x 15cm; 40cm from the top
- Absorbent Reinforcement 30cm x 40cm

Code	Description	Packaging
89543	HALYARD BASICS* Dental Drape, Sterile	Items/Case = 15




EENT FENESTRATED DRAPE, 183cm x 229cm

- 6cm x 5cm Fenestration
- CONTROL Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89421	EENT Fenestrated Drape, Sterile	Items/Case = 18

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.



MAYO STAND COVERS

Code	Description	Packaging
88665	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Sterile	Items/Case = 80
78656	Mayo Stand Cover, Reinforced HALYARD BASICS* 58cm x 137cm, Non-Sterile	Items/Case = 140
89601	Mayo Stand Cover, Reinforced, 58cm x 137cm, Sterile	Items/Case = 54
79607	Mayo Stand Cover, Reinforced, 58cm x 137cm, Non-Sterile	Items/Case = 1200
89602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
79602	Mayo Stand Cover, Reinforced, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 120
14007	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Sterile	Items/Case = 44
76358	Mayo Stand Cover, Padded, X-Large, 76cm x 145cm, Non-Sterile	Items/Case = 68
89603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Sterile	Items/Case = 44
79603	Mayo Stand Cover, Reinforced, XX-Large, 80cm x 145cm, Non-Sterile	Items/Case = 120



BACK TABLE COVERS

Code	Description	Packaging
88666	Back Table Cover, Zone Reinforced HALYARD BASICS* 112cm x 229cm, Sterile	Items/Case = 45
78664	Back Table Cover, Zone Reinforced HALYARD BASICS* 112cm x 229cm, Non-Sterile	Items/Case = 60
89632	Back Table Cover, HALYARD BASICS*, 100cm x 150cm, Sterile	Items/Case = 50
79482	Back Table Cover, 115cm x 152cm, Non-Sterile, Standard Fold	Items/Case = 90
42217	Back Table Cover, Standard, 112cm x 229cm, Sterile	Items/Case = 28
42218	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Unfolded	Items/Case = 640
42231	Back Table Cover, Standard, 112cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42224	Back Table Cover, Standard, 152cm x 229cm, Sterile	Items/Case = 36
42225	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 480
42233	Back Table Cover, Standard, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 400
42306	Back Table Cover, Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 28
42307	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 400
42366	Back Table Cover, Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 200
42309	Back Table Cover, Heavy Duty, 177cm x 279cm, Sterile	Items/Case = 15
42310	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Unfolded	Items/Case = 280
42367	Back Table Cover, Heavy Duty, 177cm x 279cm, Non-Sterile, Fan Fold	Items/Case = 160
29258	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Sterile	Items/Case = 20
29259	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Fan Fold	Items/Case = 22
29269	Back Table Cover, Padded Heavy Duty, 152cm x 229cm, Non-Sterile, Unfolded	Items/Case = 300
37498	Back Table Cover, Two Tier Heavy Duty, 152cm, Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

LAPAROSCOPY PACKS




LAPAROSCOPY PACK I

Code	Description	Packaging
88088	Laparoscopy Pack I	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm 1-Absorbent Towel, 38cm x 56cm
89411		1-Under Buttocks Drape, Plastic, 75cm x 112cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
89119	79117	1-Laparoscopy "T" Drape with CONTROL* Plus Reinforcement, 249cm x 104cm x 251cm <ul style="list-style-type: none"> • 15cm x 20cm Rectangular Abdominal Fenestration • 7cm x 15cm Perineal Fenestration with Flap • Tube Holders



LAPAROSCOPY / PELVISCOPY PACK

Code	Description	Packaging
88222	Laparoscopy/Pelviscopy Pack	Items/Case = 6
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large 1-ULTRA Surgical Gown, with Towel, Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
89413		1-Under Buttocks Drape, Plastic-Reinforced, Large 101cm x 116cm
89408	79408	2-Leggings, with 15cm Cuff, Sterile Pair, 78cm x 122cm 1-Laparoscopy / Pelviscopy Drape with CONTROL* Plus Reinforcement (without Leggings), 264cm x 89cm x 264cm <ul style="list-style-type: none"> • 18cm x 23cm Trapezoidal Abdominal Fenestration • 8cm x 15cm Perineal Fenestration with Flap



LAPAROSCOPY / PELVISCOPY DRAPE, 254cm x 76cm x 312cm

- 17cm x 23cm Trapezoidal Abdominal Fenestration
- 7cm x 20cm Perineal Fenestration with Flap
- CONTROL* Plus Fabric Reinforcement
- Attached Leggings
- Hook and Loop Tube and Line Controls

Code	Description	Packaging
89224	Laparoscopy/Pelviscopy Drape, Sterile	Items/Case = 10
79224	Laparoscopy/Pelviscopy Drape, Non-Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

LITHOTOMY PACKS



LITHOTOMY PACK II

Code	Description	Packaging
88521	Lithotomy Pack II	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
		1-Absorbent Towel, 38cm x 56cm
89413	79413	1-Under Buttocks Drape, Plastic-Reinforced, 102cm x 117cm
89408	79408	2-Leggings with 15cm Cuff, 78cm x 122cm
47614	79101	1-Half Drape, 147cm x 101cm



LITHOTOMY PACK V

Code	Description	Packaging
88551	Lithotomy Pack V	Items/Case = 12
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
89701	79720	1-Absorbent Towel, 38cm x 56cm
89408	79408	2-Leggings, with 15cm Cuffs
	79551	1-Lithotomy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 78cm x 190cm with 8cm x 15cm Perineal Fenestration

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


LITHOTOMY DRAPES



LITHOTOMY DRAPE, 254cm x 76cm x 221cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Attached Leggings


Code	Description	Packaging
89203	Lithotomy Drape, Sterile	Items/Case = 16



GYN / UROLOGY DRAPE WITH POUCH, 160cm x 284cm x 332cm

- Pouch with Mesh Screen and Port
- Attached Leggings
- CONTROL* Plus Fabric Reinforcement
- 7cm x 15cm Perineal Fenestration
- Tube Holders
- Armboard Covers


Code	Description	Packaging
89204	Gyn/Urology Drape with Pouch, Sterile	Items/Case = 10
79204	Gyn/Urology Drape with Pouch, Non-Sterile	Items/Case = 10



LITHOTOMY "T" DRAPE, 152cm x 78cm x 190cm

- 7cm x 15cm Perineal Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
79551	Lithotomy "T" Drape, Non-Sterile	Items/Case = 55



LEGGINGS, 78cm x 122cm, WITH 15cm CUFF

Code	Description	Packaging
89408	Leggings, Sterile Pair	Items/Case = 26
79408	Leggings, Non-Sterile	Items/Case = 70

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


UNDER BUTTOCKS DRAPES



UNDER BUTTOCKS DRAPE, ZONE-REINFORCED, MEDIUM, 96cm x 112cm


- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
79411	Under Buttocks Drape, Zone-Reinforced, Medium, Non-Sterile	Items/Case = 115



UNDER BUTTOCKS DRAPE, PLASTIC REINFORCED, LARGE, 101cm x 116cm


Code	Description	Packaging
89413	Under Buttocks Drape, Plastic Reinforced, Large, Sterile	Items/Case = 50



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- CONTROL* Plus Absorbent Pad
- 72cm x 66cm Triangular Pouch

Code	Description	Packaging
89414	Under Buttocks Drape, with Fluid Collection Pouch, Sterile	Items/Case = 30
79414	Under Buttocks Drape, with Fluid Collection Pouch, Non-Sterile	Items/Case = 90



UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH, 101cm x 112cm

- Screen and Suction Port
- Adhesive Pouch Closure
- 69cm x 61cm x 7cm Triangular Pouch

Code	Description	Packaging
89415	Under Buttocks Drape with Fluid Collection Pouch, Sterile	Items/Case = 40
79415	Under Buttocks Drape with Fluid Collection Pouch, Non-Sterile	Items/Case = 100

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.



MINOR PROCEDURE FENESTRATED DRAPE, 76cm x 76cm

- 4cm x 5cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89206	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 100
79206	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 250



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 7cm Diameter Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89207	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79207	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 122cm x 152cm

- 12cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape


Code	Description	Packaging
89208	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 40
79208	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 60



MINOR PROCEDURE FENESTRATED DRAPE, 193cm x 315cm

- 15cm x 15cm Fenestration
- CONTROL* Fabric Reinforcement
- Positioning Tape

Code	Description	Packaging
89209	Minor Procedure Fenestrated Drape, Sterile	Items/Case = 19
79209	Minor Procedure Fenestrated Drape, Non-Sterile	Items/Case = 23



CIRCUMCISION DRAPE, 31cm x 51cm

- 2.5cm Circular Fenestration
- CONTROL* Fabric Reinforcement

Code	Description	Packaging
44972	Circumcision Drape, Sterile	Items/Case = 40

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.



**BASIC OPHTHALMIC DRAPE,
132cm x 142cm**

- Lightweight Fabric
- 12cm x 9cm Adhesive Incise Area
- Fluid Collection Pouch, with Eye Wick

Code	Description	Packaging
89304	Basic Ophthalmic Drape, Sterile	Items/Case = 40
79304	Basic Ophthalmic Drape, Non-Sterile	Items/Case = 100




**OCUARTS* FULL BODY OPHTHALMIC DRAPE,
183cm x 292cm**

- 11cm x 11cm Adhesive Incise Area
- 5cm x 3cm Oval Fenestration
- 12cm x 15cm Optional Incise Drape
- 270° Fluid Collection Pouch, Foam Band
- Breather Bar
- Tube Holders


Code	Description	Packaging
89356	OCUARTS* Full Body Ophthalmic Drape, Sterile	Items/Case = 10
79356	OCUARTS* Full Body Ophthalmic Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.



ORTHO PACK I – MAJOR

Code	Description	Packaging
88411	Ortho Pack I – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89618	79498	1-Back Table Cover, Zone-Reinforced, 112cm x 228cm
44682	44664	2-AERO CHROME* Breathable High Performance Surgical Gown, with Towel, X-Large
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm
89331	79330	1-Suture Bag 1-Impervious Split Drape, 152cm x 178cm with 10cm x 53cm Split
89381	79360	1-Ortho U-Drape with CONTROL* Plus Reinforcement, 223cm x 274cm with 10cm x 102cm Split 1-Ortho Bar Drape with CONTROL* Plus Reinforcement, 254cm x 152cm



ORTHOARTS* LOWER EXTREMITY PACK

Code	Description	Packaging
88471	ORTHOARTS* Lower Extremity Pack	Items/Case = 3
Sterile Component Code	Non-Sterile Component Code	Contents
89906	79006	1-Back Table Cover, Heavy-Duty, 152cm x 228cm
95111	71311	1-ULTRA Surgical Gown, In Overwrap, with Towel, Large 1-Extremity Prep Pad
95221	74921	1-ULTRA Fabric-Reinforced Surgical Gown, with Towel, X-Large
44686	44668	1-AERO CHROME* Breathable High Performance Specialty Surgical Gown, X-Large, X-Long
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag 1-304cm Suction Tubing
89821	79515	1-Stockinette, Large
89331	79330	1-Split Drape, 152cm x 178cm with 10cm x 53cm Split
89278	13233	1-ORTHOARTS* Lower Extremity Drape with CONTROL* Plus Reinforcement, with 6cm Conformable Circular Fenestration
89053		3-Drape Tape Strips, Sterile

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


EXTREMITY DRAPE



LOWER EXTREMITY DRAPE, 274cm x 203cm x 315cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89016	Lower Extremity Drape, Sterile	Items/Case = 12



LOWER EXTREMITY DRAPE, 223cm x 332cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
89276	Lower Extremity Drape, Sterile	Items/Case = 14



ORTHOARTS* LOWER EXTREMITY DRAPE, 289cm x 223cm x 320cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89278	ORTHOARTS* Lower Extremity Drape, Sterile	Items/Case = 7
13233	ORTHOARTS* Lower Extremity Drape, Non-Sterile	Items/Case = 84



ORTHOARTS* EXTRA LONG LOWER EXTREMITY DRAPE, 289cm x 223cm x 381cm

- Armboard Covers
- 6cm Conformable Circular Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
29219	ORTHOARTS* Extra Long Lower Extremity Drape, Sterile	Items/Case = 8
29220	ORTHOARTS* Extra Long Lower Extremity Drape, Non-Sterile	Items/Case = 8

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


EXTREMITY DRAPES



UNIVERSAL EXTREMITY DRAPE, 229cm x 333cm

- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration


Code	Description	Packaging
89279	Universal Extremity Drape, Sterile	Items/Case = 14
79279	Universal Extremity Drape, Non-Sterile	Items/Case = 150



BILATERAL LIMB DRAPE, 193cm x 305cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders


Code	Description	Packaging
89291	Bilateral Limb Drape, Sterile	Items/Case = 14



BILATERAL LIMB DRAPE WITH ARMBOARD COVERS, 295cm x 223cm x 366cm

- CONTROL* Plus Fabric Reinforcement
- Two 6cm Conformable Circular Fenestrations
- Tube Holders
- Armboard Covers

Code	Description	Packaging
89380	Bilateral Limb Drape with Armboard Covers, Sterile	Items/Case = 12



UPPER EXTREMITY DRAPE, 193cm x 305cm

- 3cm Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89271	Upper Extremity Drape, Sterile	Items/Case = 14
79273	Upper Extremity Drape, Non-Sterile	Items/Case = 20

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


HAND DRAPES



HAND DRAPE,
289cm x 360cm

- 157cm x 147cm Hand Table Cover
- 3cm Conformable Circular Fenestration
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
77450	Hand Drape, Sterile	Items/Case = 11
13227	Hand Drape, Non-Sterile	Items/Case = 108




ORTHOARTS* HAND DRAPE,
193CM X 381CM

- 157cm x 147cm Gusseted Hand Table Cover
- 3cm Conformable Fenestration
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89272	ORTHOARTS* Hand Drape, Sterile	Items/Case = 6

All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.


KNEE ARTHROSCOPY



KNEE ARTHROSCOPY DRAPE, 228cm x 315cm

- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 8cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders


Code	Description	Packaging
89265	Knee Arthroscopy Drape, Sterile	Items/Case = 10
79265	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 150



KNEE ARTHROSCOPY DRAPE WITH ARMBOARD COVERS, 290cm x 223cm x 315cm

- Armboard Covers
- Fluid Collection Pouch with Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89371	Knee Arthroscopy Drape with Armboard Covers, Sterile	Items/Case = 10
79371	Knee Arthroscopy Drape with Armboard Covers, Non-Sterile	Items/Case = 13



ORTHOARTS* KNEE ARTHROSCOPY DRAPE, 289cm x 223cm x 294cm

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 76cm x 99cm x 9cm Triangular Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
89267	ORTHOARTS* Knee Arthroscopy Drape, Sterile	Items/Case = 7
79267	ORTHOARTS* Knee Arthroscopy Drape, Non-Sterile	Items/Case = 7

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.


These products are not available for purchase by the general public.

KNEE ARTHROSCOPY



HALYARD BASICS* KNEE ARTHROSCOPY PACK


Code	Description	Packaging
88886	HALYARD BASICS* Knee Arthroscopy Pack	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
88660	78660	1-HALYARD BASICS* Film Split Drape, 152cm x 178cm, with 10cm x 53cm Split
88661	78661	1-HALYARD BASICS* Knee Arthroscopy Drape with Fluid Collection Pouch, 229cm x 315cm
89564	79764	1-HALYARD BASICS* Back Table Cover, Zone-Reinforced, 152cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm



ORTHOARTS* EXTRA LONG KNEE ARTHROSCOPY DRAPE, 294cm x 233cm x 447cm

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- 95cm x 112cm Pentagonal Pouch
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
36930	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Sterile	Items/Case = 10
42429	ORTHOARTS* Extra Long Knee Arthroscopy Drape, Non-Sterile	Items/Case = 10



HALYARD BASICS* KNEE ARTHROSCOPY DRAPE, 229cm x 315cm

- Fluid Collection Pouch with Suction Port
- 6cm Conformable Circular Fenestration
- Tube Holders

Code	Description	Packaging
88661	Knee Arthroscopy Drape, Sterile	Items/Case = 20
78661	Knee Arthroscopy Drape, Non-Sterile	Items/Case = 30

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


SHOULDER DRAPES



ORTHOARTS* BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE, 411cm x 261cm

- Attached Fluid Collection Pouch
- 95cm x 97cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Heavy Duty Fabric


Code	Description	Packaging
36932	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
42426	ORTHOARTS* Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 10



BEACH CHAIR SHOULDER DRAPE WITH ATTACHABLE FLUID COLLECTION POUCH, 437cm x 254cm

- Attachable Fluid Collection Pouch with Suction Port, 89cm x 109cm
- CONTROL* Plus Fabric Reinforcement
- Hook and Loop Tube and Line Control
- Pouch Attachment Zones
- 90cm x 110cm Trapezoidal Pouch

Code	Description	Packaging
89066	Beach Chair Shoulder Drape with Attachable Fluid Collection Pouch, Sterile	Items/Case = 6
79049	Beach Chair Shoulder Drape without Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 11
79061	Attachable Fluid Collection Pouch, Non-Sterile	Items/Case = 44



BEACH CHAIR SHOULDER ARTHROSCOPY DRAPE II, 406cm x 259cm

- Attached Fluid Collection Pouch
- 103m x 63cm x 2cm Pentagonal Pouch
- Hook and Loop Tube and Line Control
- 6cm Conformable Tear Drop Fenestration
- Instrument Pad Fabric Reinforcement

Code	Description	Packaging
89070	Beach Chair Shoulder Arthroscopy Drape, Sterile	Items/Case = 10
79070	Beach Chair Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric



Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

SHOULDER DRAPES



2-PIECE SHOULDER ARTHROSCOPY DRAPE PACK


- Body Split Drape, 223cm x 294cm
- 10cm x 68cm Split
- CONTROL* Plus Fabric Reinforcement
- Shoulder Arthroscopy Drape, 254cm x 157cm
- Fluid Collection Pouch with Suction Ports
- 55cm x 112cm Trapezoidal Pouch
- Tube Holders

Code	Description	Packaging
89490	2-Piece Shoulder Arthroscopy Drape Pack, Sterile	Items/Case = 12
79490	2-Piece Shoulder Arthroscopy Drape Pack, Non-Sterile	Items/Case = 12

ORTHOARTS* SHOULDER ARTHROSCOPY DRAPE, 254cm x 152cm

- Fluid Collection Pouch
- 79cm x 149cm x 48cm Trapezoidal Pouch
- Heavy-Duty Fabric


Code	Description	Packaging
89492	ORTHOARTS* Shoulder Arthroscopy Drape, Sterile	Items/Case = 12
79492	ORTHOARTS* Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 15



ORTHOARTS* SHOULDER SPLIT DRAPE, 284cm x 152cm

- 15cm x 20cm Dual Splits
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders


Code	Description	Packaging
89274	ORTHOARTS* Shoulder Split Drape, Sterile	Items/Case = 12



HALYARD BASICS* SHOULDER SPLIT DRAPE, 285cm x 152cm

- 13cm x 34cm Split
- Fabric Reinforcement


Code	Description	Packaging
89559	Shoulder Split Drape, Sterile	Items/Case = 30
79759	Shoulder Split Drape, Non-Sterile	Items/Case = 50



HALYARD BASICS* SHOULDER ARTHROSCOPY DRAPE, 254cm x 158cm

- Fluid Collection Pouch with Suction Port
- Tube Holders

Code	Description	Packaging
89558	Shoulder Arthroscopy Drape, Sterile	Items/Case = 30
79758	Shoulder Arthroscopy Drape, Non-Sterile	Items/Case = 50



All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


U & BAR DRAPES



U-DRAPE, 193cm x 305cm

- 10cm x 101cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
89301	U-Drape, Sterile	Items/Case = 17
13231	U-Drape, Non-Sterile	Items/Case = 190



LARGE U-DRAPE, 223cm x 317cm

- 15cm x 137cm Split
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape


Code	Description	Packaging
89311	Large Ortho U-Drape, Sterile	Items/Case = 11
79363	Large Ortho U-Drape, Non-Sterile	Items/Case = 15



ORTHOARTS* U-DRAPE, 223CM X 274CM

- 10cm x 101cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement

Code	Description	Packaging
89347	ORTHOARTS* U-Drape, Sterile,	Items/Case = 10
13294	ORTHOARTS* U-Drape, Non-Sterile	Items/Case = 132



ORTHOARTS* BODY SPLIT DRAPE, 223CM X 294CM

- 10cm x 68cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
89275	ORTHOARTS* Body Split Drape, Sterile	Items/Case = 10
79275	ORTHOARTS* Body Split Drape, Non-Sterile	Items/Case = 12

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


U & BAR DRAPES



SPLIT DRAPE, 152cm x 178cm

- Polyethylene Material
- Impervious
- 10cm x 53cm Split


Code	Description	Packaging
89331	Split Drape, Sterile	Items/Case = 42
79330	Split Drape, Non-Sterile	Items/Case = 115



LARGE SPLIT DRAPE, 152cm x 205cm

- Polyethylene Material
- Impervious
- 11cm x 78cm Split


Code	Description	Packaging
89333	Large Split Drape, Sterile	Items/Case = 34
79333	Large Split Drape, Non-Sterile	Items/Case = 100



HALYARD BASICS* BODY SPLIT DRAPE, 224cm x 295cm

- 10cm x 67cm Split
- Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89557	Body Split Drape, Sterile	Items/Case = 25
79757	Body Split Drape, Non-Sterile	Items/Case = 30



ORTHO SPLIT "T" DRAPE, 297cm x 210cm x 342cm

- Armboard Covers
- 10cm x 53cm Split
- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape

Code	Description	Packaging
36931	Ortho Split "T" Drape, Sterile	Items/Case = 10

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric


Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.


U & BAR DRAPES



ORTHO BAR DRAPE, 254cm x 152cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders

Code	Description	Packaging
89381	Ortho Bar Drape, Sterile	Items/Case = 18
79360	Ortho Bar Drape, Non-Sterile	Items/Case = 28



VERTICAL ISOLATION DRAPE, 228cm x 320cm

- Clear Film
- Adhesive Incise Area 40cm x 76cm
- Attachable Fluid Collection Pouch with Suction Port
- 42cm x 84cm x 17cm Trapezoidal Pouch
- Tube Holders
- Central Release Tape Strip For Horizontal Rod Attachment

Code	Description	Packaging
77723	Vertical Isolation Drape, Sterile	Items/Case = 14

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.




THYROID DRAPE,
193cm x 315cm

- 12cm x 12cm Diamond Fenestration
- CONTROL* Plus Fabric Reinforcement
- Instrument Pad
- Tube Holders

Code	Description	Packaging
89255	Thyroid Drape, Sterile	Items/Case = 13
79255	Thyroid Drape, Non-Sterile	Items/Case = 18


All codes in **bold** indicate a sterile product.
† Reinforcements include:
CONTROL* Plus – Maximum Absorbency
CONTROL* – Absorbent Reinforcement Fabric
Instrument Pad – Non-Slip Pad
NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.
These products are not available for purchase by the general public.



TOP DRAPE, 273cm x 127cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759


Code	Description	Packaging
34906	Top Drape, Sterile	Items/Case = 15
34901	Top Drape, Non-Sterile	Items/Case = 275



BOTTOM DRAPE, 193cm x 195cm

- CONTROL* Plus Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34902	Bottom Drape, Sterile	Items/Case = 12
34903	Bottom Drape, Non-Sterile	Items/Case = 242



SIDE DRAPE, 195cm x 105cm

- CONTROL* Plus Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88759

Code	Description	Packaging
34904	Side Drape, Sterile	Items/Case = 20
34905	Side Drape, Non-Sterile	Items/Case = 440

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.



TOP DRAPE, 276cm x 150cm


- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769

Code Description

79095 Top Drape, Non-Sterile

Packaging

Items/Case = 60



BOTTOM DRAPE, 165cm x 193cm


- CONTROL* Fabric Reinforcement
- Tube Holders
- Positioning Tape
- From Universal Pack 88769

Code Description

79096 Bottom Drape, Non-Sterile

Packaging

Items/Case = 70



SIDE DRAPE, 114cm x 84cm

- CONTROL* Fabric Reinforcement
- Positioning Tape
- From Universal Pack 88769

Code Description

79097 Side Drape, Non-Sterile

Packaging

Items/Case = 200

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.



UNIVERSAL PACK – MAJOR WITH LONG SIDE DRAPES

Code	Description	Packaging
88759	Universal Pack – Major	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Fully Reinforced, 137cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Absorbent Towel, 15cm x 56cm 1-Suture Bag
34902	34903	1-Bottom Drape with CONTROL* Plus Reinforcement, 193cm x 195cm
34906	34901	1-Top Drape with CONTROL* Plus Reinforcement, 273cm x 127cm
34904	34905	2-Side Drapes with CONTROL* Plus Reinforcement, 195cm x 105cm



UNIVERSAL PACK – MINOR WITH SHORT SIDE DRAPES

Code	Description	Packaging
88769	Universal Pack – Minor	Items/Case = 14
Sterile Component Code	Non-Sterile Component Code	Contents
89611	79613	1-Back Table Cover, Fully-Reinforced, 112cm x 224cm
89601	79601	1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag
	79097	2-Side Drapes with CONTROL* Reinforcement, 114cm x 84cm
	79096	1-Bottom Drape with CONTROL* Reinforcement, 165cm x 193cm
	79095	1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.



UNIVERSAL PACK I

Code	Description	Packaging
88878	Universal Pack – I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
88666	78658	1-HALYARD BASICS* Back Table Cover, Zone Reinforced, 112cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm
	79756	2-HALYARD BASICS* Side Drape, 193 x 112cm
	79754	1-HALYARD BASICS* Top Drape, 274 x 127cm
	79755	1-HALYARD BASICS* Bottom Drape, 193 x 173cm



UNIVERSAL PACK III WITH MEDIUM SIDE DRAPES

Code	Description	Packaging
88774	Universal Pack – III	Items/Case = 8
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Fully-Reinforced, 133cm x 248cm 1-Mayo Stand Cover, 58cm x 137cm 1-Suture Bag 2-Side Drapes with CONTROL* Reinforcement, 145cm x 99cm 1-Bottom Drape with CONTROL* Reinforcement, 193cm x 193cm 1-Top Drape with CONTROL* Reinforcement, 276cm x 150cm 2-Absorbent towels, 38cm x 56cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.



U-BAR* PACK

Code	Description	Packaging
88334	U-BAR* Pack	Items/Case = 16
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Zone Reinforced, 112cm x 224cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-Suture Bag 1-U-Drape with CONTROL* Reinforcement, 180cm x 315cm, with 10cm x 102cm Split 1-Bar Drape with CONTROL* Reinforcement, 105cm x 228cm



U-BAR* PACK – IV

Code	Description	Packaging
88773	U-BAR* Pack – IV	Items/Case = 10
Sterile Component Code	Non-Sterile Component Code	Contents
89601	79601	1-Back Table Cover, Reinforced, 133cm x 248cm 1-Mayo Stand Cover, Reinforced, 58cm x 137cm 1-U-Drape with CONTROL* Reinforcement, 192cm x 315cm, with 11cm x 51cm Split 1-Bar Drape with CONTROL* Reinforcement, 248cm x 145cm

HALYARD BASICS* U-BAR PACK I

Code	Description	Packaging
88875	HALYARD BASICS* U-BAR Pack I	Items/Case = 5
Sterile Component Code	Non-Sterile Component Code	Contents
89532	79532	1-HALYARD BASICS* U-Drape, 203cm x 315cm, with 10cm x 112cm Split
89533	79533	1-HALYARD BASICS* Bar Drape, 244cm x 112cm
88666	78658	1-HALYARD BASICS* Back Table Cover, Zone-Reinforced, 112cm x 229cm
88665	78656	1-HALYARD BASICS* Mayo Stand Cover, Reinforced, 58cm x 137cm

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency


CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.



CYSTOSCOPY PACK I

Code	Description	Packaging
88611	Cystoscopy Pack I	Items/Case = 12
Sterile Component Code	Non-Sterile Component Code	Contents
89408	79408	1-Back Table Cover, Zone-Reinforced, 112cm x 198cm 1-Absorbent Towel, 38cm x 56cm 2-Leggings, with 15cm Cuff, 78cm x 122cm 1-Cystoscopy "T" Drape with CONTROL* Plus Reinforcement, 152cm x 76cm x 213cm with 7cm x 15cm Perineal Fenestration and 23cm x 23cm Screen

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.


Always read the label and follow the directions for use.

These products are not available for purchase by the general public.




ABSORBENT TOWELS, 38cm x 56cm

Code	Description	Packaging
89701	Absorbent Towels, White, Sterile, Pack of 2	Items/Case = 100
79720	Absorbent Towels, White, Non-Sterile	Items/Case = 1000



TAPE STRIPS, 10cm x 55cm


Code	Description	Packaging
89053	Tape Strips, Sterile, Pack of 2	Items/Case = 300



DIATHERMY POUCH, 38cm x 41cm

- Adhesive Backing Strip

Code	Description	Packaging
89945	Diathermy Pouch, Sterile	Items/Case = 150
79945	Diathermy Pouch, Non-Sterile	Items/Case = 220



UTILITY DRAPE WITH TAPE, 38cm x 64cm

Code	Description	Packaging
89731	Utility Drape with Tape, Sterile, Pack of 2	Items/Case = 100
79731	Utility Drape with Tape, Non-Sterile	Items/Case = 550



STOCKINETTES

Code	Description	Packaging
89811	Stockinette, Sterile, Medium, 23cm x 112cm	Items/Case = 36
79505	Stockinette, Non-Sterile, Medium, 24cm x 117cm	Items/Case = 75
89821	Stockinette, Sterile, Large, 30cm x 122cm	Items/Case = 25
79515	Stockinette, Non-Sterile, Large, 29cm x 117cm	Items/Case = 50

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad

NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.

These products are not available for purchase by the general public.

O.R. ACCESSORIES



HALYARD* STIRRUP STRAP

- Designed to ensure patient comfort in the lithotomy position
- Made of soft fabric and featuring a layer of padding
- Comes with slip ring for easy adjustment
- Single-patient use

Code	Description	Size	Packaging
52712	Personal Stirrup Strap with Slip Ring	48cm x 9cm	10 eaches/pack 4 packs/case 40 eaches/case



HALYARD* TS-30* PERSONAL UTILITY STRAP

- Designed to secure and position the patient during a surgical procedure
- Hook and loop closure for easy adjustment of the strap
- Soft fabric for patient comfort
- Single-patient use

Code	Description	Size	Packaging
55299	TS-30* Personal Utility Strap	76cm x 8cm	12 eaches/pack 2 packs/case 24 eaches/case



HALYARD* TOURNIQUET COVER

- Designed to cover a reusable, pneumatic tourniquet cuff to protect it from skin prep and other procedural fluids
- Soft padding helps prevent pinching and wrinkling of patients' skin
- Comes in two sizes – small and large
- Single-patient use

Code	Description	Size	Packaging
53600	O.R. Tourniquet Cover, Small	53.5cm x 37cm	15 eaches/pack 2 packs/case 30 eaches/case
53605	O.R. Tourniquet Cover, Large	78cm x 37cm	15 eaches/pack 2 packs/case 30 eaches/case

All codes in **bold** indicate a sterile product.

† Reinforcements include:

CONTROL* Plus – Maximum Absorbency

CONTROL* – Absorbent Reinforcement Fabric

Instrument Pad – Non-Slip Pad


NOTE: Please see package inserts for directions for use for each of these products.

Always read the label and follow the directions for use.



DAY SURGERY PRODUCTS

SURGICAL GOWNS



SURGICAL GOWNS

THE ASSOCIATION FOR THE ADVANCEMENT OF MEDICAL INSTRUMENTATION (AAMI)

DO YOUR GOWNS COMPLY?



THE AAMI STANDARD: AAMI PB70:2012

Liquid Barrier Performance and Classification of Protective Apparel and Drapes Intended for Use in Health Care Facilities

The AAMI Standard classifies the barrier performance of surgical gowns through the use of four industry standard tests:

AAMI:PB70	Risk of Exposure: Fluid Amount	Risk of Exposure: Fluid Spray or Splash	Risk of Exposure: Pressure on Gown
Level 1	Minimal	Minimal	Minimal
Level 2	Low	Low	Low
Level 3	Moderate	Moderate	Moderate
Level 4	High	High	High

LEVEL	TEST†
1	Impact penetration test (Gowns & Drapes)
2	Impact penetration + Hydrostatic pressure test = $\geq 20\text{cmH}_2\text{O}$ (Gowns & Drapes)
3	Impact penetration + Hydrostatic pressure test = $\geq 50\text{cmH}_2\text{O}$ (Gowns & Drapes)
4	ASTM F1670 – Standard test method for Resistance of Materials used in protective clothing to Penetration by synthetic blood (Drapes only)
	ASTM F1671 – Standard test method for Resistance of Materials used in protective clothing to Penetration by Blood-Borne Pathogens (Gowns only)

† Applies to products' critical zones

SURGICAL GOWNS

ASTM

(AMERICAN SOCIETY FOR TESTING AND MATERIALS)

LEVEL 4



ASTM F1670 – TESTS FOR FLUID PENETRATION

Standard Test Method for Resistance of Materials used in Protective Clothing to Penetration by Synthetic Blood

- If it passes this test it is considered impervious to blood

ASTM F1671 – TESTS FOR VIRAL PENETRATION

Standard Test Method for Resistance of Materials used in Protective Clothing to Penetration by Blood-borne Pathogens using Phi-X174 Bacteriophage Penetration as a Test System

- Pass/Fail Tests
- Essential to meet level 4 under AAMI PB70:2012

HALYARD GOWNS

HALYARD* AERO CHROME* BREATHABLE PERFORMANCE AAMI 4 SURGICAL GOWNS

**What makes AERO CHROME*
an AAMI 4 gown?**

The critical zones of the gown consist of a fabric that is highly effective in preventing:

1. Blood penetration
2. Penetration of blood-borne viruses
3. Splash penetration
4. Fluid penetration when exposed to high pressures ($\geq 50\text{cm H}_2\text{O}$).

**AERO CHROME* is ideal
for high level of fluid and
microbial protection in
the critical zones.**



HALYARD* AERO BLUE* PERFORMANCE AAMI 3 SURGICAL GOWNS

**What makes AERO BLUE*
an AAMI 3 gown?**

The critical zones of the gown consist of a fabric that is highly effective in preventing:

1. Splash penetration
2. Fluid penetration when exposed to high pressures (up to $400\text{cm H}_2\text{O}$ ²² versus the standard of $50\text{cm H}_2\text{O}$).

**AERO BLUE* is ideal for moderate
fluid procedures.**

NOT ALL AAMI LEVEL 3 GOWNS ARE CREATED EQUAL

AERO BLUE* offers the highest protection available in an AAMI Level 3 gown:²³

**Four times more fluid protection in the
critical zones than other gowns in its class.¹⁹**





AERO CHROME* BREATHABLE PERFORMANCE SURGICAL GOWN

- Soft Knitted V-neck Collar
- Clearly Identifiable Silver Colour
- Soft, Light and Breathable Fabric
- Raglan Sleeves
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes
- AAMI Level 4 / EN 13795 High Performance


Code	Description	Packaging
44680	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Small	Items/Case = 30
44681	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large	Items/Case = 30
44682	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
44683	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
44684	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24
44685	AERO CHROME* Breathable Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
44686	AERO CHROME* Breathable Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
44687	AERO CHROME* Breathable Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28



AERO BLUE* PERFORMANCE SURGICAL GOWNS

- Patented Breathable Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour-Coded Neck Binding – Green
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes
- AAMI Level 3 / EN 13795 High Performance

Code	Description	Packaging
41870	AERO BLUE* Performance Surgical Gown, Sterile, Small	Items/Case = 30
41871	AERO BLUE* Performance Surgical Gown, Sterile, Large	Items/Case = 30
42620	AERO BLUE* Performance Surgical Gown, Sterile, Large, X-Long	Items/Case = 28
41872	AERO BLUE* Performance Surgical Gown, Sterile, X-Large	Items/Case = 28
42621	AERO BLUE* Performance Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 30
41873	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large	Items/Case = 26
43940	AERO BLUE* Performance Surgical Gown, Sterile, XX-Large, X-Long	Items/Case = 28
41874	AERO BLUE* Performance Surgical Gown, Sterile, XXX-Large	Items/Case = 24



SURGICAL GOWNS


- EVOLUTION* 4 Fabric
- Low Linting Fabric
- Set-In-Sleeves
- Colour-Coded Neck Binding – Yellow
- Adjustable Hook and Loop Neckline
- Book Fold
- 1 x Overwrap, 90cm x 90cm
- 1 x Blue Huck Towel Included for all Sterile Codes, 51cm x 64cm
- EN 13795 Standard Performance

Code	Description	Packaging
90003	Surgical Gown, Sterile, Small to Medium	Items/Case = 36
90018	Surgical Gown, Sterile, Large	Items/Case = 36
90019	Surgical Gown, Sterile, Large, X-Long	Items/Case = 34
90048	Surgical Gown, Sterile, X-Large	Items/Case = 32
90049	Surgical Gown, Sterile, X-Large, X-Long	Items/Case = 34
90050	Surgical Gown, Sterile, XX-Large	Items/Case = 28
90138	Surgical Gowns Double Pack	Items/Case = 18
	Contents:	
72978	1-Surgical Gown, Large	
72988	1-Surgical Gown, X-Large	



DAY SURGERY PRODUCTS

CUSTOM PROCEDURE PACKS (CPP)



CUSTOM PROCEDURE PACKS (CPP)

THE FUTURE OF CUSTOM PROCEDURE PACKS FROM THE NAME YOU TRUST

For over 30 years HALYARD has been providing high quality surgical drapes and gowns, designed to deliver protection, comfort and optimal patient outcomes.

Leveraging Halyard's global expertise in supply chain, component sourcing, utilisation reviews and cost efficiencies, we are now extending our offering to provide ANZ clinicians with HALYARD branded Custom Procedure Packs (CPPs).

Based on years of working together with clinicians globally we know that some of the challenges you face in Operating Theatres include:

- Operational inefficiencies
- Lack of visibility to costs and inventories
- Suboptimal utilisation of Operating Theatre supplies
- Procedural supply waste
- Need for improved throughput
- Staff satisfaction and retention


That's where HALYARD comes in. When you make the decision to switch to HALYARD CPPs, our clinical experts work with you to ensure EVERY component is used EVERY time. It's part of our PACT pledge and the HALYARD CPP Solution.





DAY SURGERY PRODUCTS

PERSONAL PROTECTIVE EQUIPMENT (PPE)



PERSONAL PROTECTIVE EQUIPMENT

Halyard offers products across all areas of Personal Protective Equipment (PPE) where there may be risk of exposure to blood or other bodily fluids.

SEQUENCE FOR **DONNING** PPE:

1. PROTECT YOUR BODY

Ensure the gown:

- Covers torso from neck to knees
- Covers arms to end of wrists
- Wraps around the sides to fasten at neck and back



2. PROTECT YOUR FACE

- Create a curve in the flexible band for the nose bridge and expand the mask
- Secure the ties or elastics at the middle of the head and neck, or fasten the ear loops
- Fit it snugly to your face and below your chin – fit check the N95 and P2 mask
- Place the face shield or eyewear over the face and eyes and adjust to fit



3. PROTECT YOUR HANDS

- Extend the cuffs to cover the wrist of the gown
- Limit the surfaces touched
- Change gloves when torn or heavily contaminated
- Perform hand hygiene as required



FINALLY **STOP** and **CHECK** all PPE is in place



SEQUENCE FOR **ADVANCED DOFFING** PPE:

Remember! The outside and front of your PPE, including the sleeves of the gown, are contaminated. If hands become visibly contaminated during PPE removal, wash hands before continuing.

REMOVE GOWN THEN GLOVES

Outside of gown and gloves are **CONTAMINATED!**

- 1a.** ● With gloved hands, grasp the front of the gown at the shoulders and pull the gown forward to snap the tabs at the back
- Peel off the gown, rolling it inside out
- 1b.** ● Continue over the gloves so that they too are rolled inside out
- Discard gown and gloves into a waste container
- Perform hand hygiene using ABHR or soap and water



2. REMOVE PROTECTIVE EYEWEAR

- Grasp the ear or head piece and lift away from face
- Place reusable eyewear in the designated receptacle for cleaning or
- Discard disposable eyewear into waste container for disposal



3. REMOVE MASK OR PARTICULATE RESPIRATOR (N95 OR P2)

- Remove by touching ties or elastics only, starting with the bottom then top ties or bands
- Discard in designated waste container



4. PERFORM HAND HYGIENE

- Wash hands with soap and water or use an alcohol-based hand rub immediately after removing PPE





HALYARD* ISOLATION GOWNS

- For use when expected risk of exposure to fluid is between low and moderate
- Medium-weight tri-layer fabric
- Tape-tab neck closure
- Tie waist, elastic cuffs

Code	Description	Packaging
69979	Isolation Gown, Yellow, Large	Eaches/Bundle = 10 Bundles/Case = 10 Eaches/Case = 100
69988	Isolation Gown, Yellow, X-Large	Eaches/Bundle = 10 Bundles/Case = 10 Eaches/Case = 100



HALYARD* IMPERVIOUS COMFORT GOWNS

- Potential for moderate to heavy fluid contact in front
- Tri-layer spunbond with polyethylene laminate coating
- Breathable spunbond back
- Knit cuffs
- Velcro neck tabs
- Tie waist

Code	Description	Packaging
69455	Impervious Open Back Gown, Knit Cuff, Blue, Large	Eaches/Bundle = 30 Bundles/Case = 2 Eaches/Case = 60
5511	Impervious Full Back Comfort Gown, Yellow, Regular	Eaches/Case = 50
5512	Impervious Full Back Comfort Gown, Yellow, X-Large	Eaches/Case = 50





HALYARD* IMPERVIOUS THUMBS-UP* FILM GOWNS

- Potential for soaking and heavy fluid contact in front
- Impervious polyethylene film
- Open back, thumbhooks

Code	Description	Packaging
7000	Impervious Thumbs-Up* Film Gown with Thumbhooks, Blue, Regular	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75
7001	Impervious Thumbs-Up* Film Gown with Thumbhooks, Blue, X-Large	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75
7002	Impervious Thumbs-Up* Film Gown with Thumbhooks, Yellow, Regular	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75
7003	Impervious Thumbs-Up* Film Gown with Thumbhooks, Yellow, X-Large	Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75



HALYARD* CHEMOTHERAPY PROCEDURE GOWN

- Potential for soaking and heavy fluid contact in front
- Tested against 52 chemotherapy drugs
- Meets Liquid Barrier Performance Standards ASTM 1670 and 1671
- Highest level of protection against flammability
- Polypropylene with polyethylene coating



Code	Description	Packaging
69606	Comfort Gown, Knit Cuffs, Perforated Open Back, Blue, Universal	Eaches/Case = 100
37284	Comfort Gown, Knit Cuffs, Perforated Open Back, Blue, XX-Large	Eaches/Case = 100



HALYARD* CHEMO360* PROCEDURE GOWN

- Closed in the back
- Secure tape closure

Code	Description	Packaging
47346	CHEMO360* Gown, Closed Back, Poly-Coated, Tape Closure, Blue, Universal	Eaches/Case = 100
47347	CHEMO360* Gown, Closed Back, Poly-Coated, Tape Closure, Blue, XX-Large	Eaches/Case = 100

HALYARD, YOUR PARTNER IN COMPLIANCE TO THE MASK STANDARD AS 4381:2015 FOR SINGLE-USE FACE MASKS


NEW STANDARD AS 4381:2015 (RED INDICATES CHANGE)

AS 4381: 2015 SINGLE-USE FACE MASKS

CHARACTERISTICS	LEVEL 1 1	LEVEL 2 2	LEVEL 3 3	TEST METHOD
	Level 1 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the minimum velocity specified in Table 2, bacterial filtration efficiency and differential pressure. APPLICATIONS: For general purpose medical procedures, where the wearer is not at risk of blood or bodily fluid splash or to protect staff and/or the patient from droplet exposure to microorganisms (e.g. patient with upper respiratory tract infection visits GP)	Level 2 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the middle velocity specified in Table 2, bacterial filtration efficiency and differential pressure. APPLICATIONS: For use in emergency departments, dentistry, changing dressings on small or healing wounds where minimal blood droplet exposure may possibly occur (e.g. endoscopy procedures)	Level 3 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the maximum velocity specified in Table 2, bacterial filtration efficiency and differential pressure. APPLICATIONS: For all surgical procedures, major trauma first aid or in any area where the health care worker is at risk of blood or bodily fluid splash (e.g. orthopaedic, cardiovascular procedures)	
Bacterial Filtration Efficiency (BFE) %	≥ 95%	≥ 98%	≥ 98%	ASTM F2101-14 or EN 14683:2014
Particulate Filtration Efficiency (PFE) % (0.1 µm)	Not Required	Not Required	Not Required	N/A
Differential Pressure (Delta P) mm H ₂ O/cm ²	< 4.0	< 5.0	< 5.0	EN 14683:2014
Resistance to penetration by synthetic blood (fluid resistance) min pressure in mm Hg for pass result	80mm Hg	120mm Hg	160mm Hg	ASTM F1862 / F1862M-13 or ISO 22609

For the highest levels of protection, make sure you choose HALYARD* Face Masks and Eyewear.

LEVEL 3



46727

NIOSH APPROVED N95 RESPIRATOR FACE MASKS


- Provides fluid resistance at 160mm Hg
- Tested against 7 chemotherapy drugs
- The highest level of protection available
- Recommended for Laser Procedures



Code	Description	Packaging
46727	FLUIDSHIELD* Level 3 N95 Particulate Filter Respirator and Surgical Mask, Pouch-Style with Headbands, Regular Size, Orange	Eaches/Box = 35 Boxes/Case = 6 Eaches/Case = 210
46827	FLUIDSHIELD* Level 3 N95 Particulate Filter Respirator and Surgical Mask, Pouch-Style with Headbands, Small Size, Orange	Eaches/Box = 35 Boxes/Case = 6 Eaches/Case = 210

NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.

LEVEL 3



HALYARD FLUIDSHIELD* FACE MASKS

- Provides fluid resistance at 160mm Hg
- Level 3 printed on mask nose piece and box
- Four layers of fluid resistance and bacterial challenge
- Superior comfort without compromising your protection with SO SOFT* inner lining
- Setting a new benchmark with POWERGUARD* protection



28797

Code	Description	Packaging
28797	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Ear Loops	Eaches/Box = 40 Boxes/Case = 10 Eaches/Case = 400



28800

28800	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Wraparound Visor and SO SOFT* Lining and SO SOFT* Ear Loops	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
-------	---	--



28804

28804	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with Wraparound Visor and SO SOFT* Lining	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
-------	---	--



47107

47107	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask, Pleat-Style with SO SOFT* Ear Loops, Foam Band, Orange	Eaches/Box = 40 Boxes/Case = 10 Eaches/Case = 400
-------	--	---



47147

47147	FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with SO SOFT* Ear Loops, Foam Band, Orange	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
-------	---	--



48207

48207	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask, Pleat-Style with Ties, Orange	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
-------	--	--



48247

48247	FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Ties, Orange	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
-------	---	--



48297

48297	FLUIDSHIELD* Level 3 Care Bear Surgical Mask, Pleat-Style with Ties, Pink and Blue Teddy Bear Print	Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300
-------	---	--

NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.

LEVEL 2



HALYARD FLUIDSHIELD* FACE MASKS

- Provides fluid resistance at 120mm Hg
- Level 2 printed on mask nose piece and box
- Four layers of fluid resistance and bacterial challenge
- Superior comfort without compromising your protection with SO SOFT* inner lining
- Setting a new benchmark with POWERGUARD* protection



62114

Code	Description	Packaging
62114	FLUIDSHIELD* Level 2 Fog-Free Surgical Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Ties, Blue with Orange Diamond Print	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100
62115	FLUIDSHIELD* Level 2 Fog-Free Procedure Mask with SO SOFT Ear Loops, Blue	Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100



62115

LEVEL 1



HALYARD FLUIDSHIELD* FACE MASKS

- Provides fluid resistance at 80mm Hg
- Level 1 printed on mask nose piece and box
- Four layers of fluid resistance and bacterial challenge
- Superior comfort without compromising your protection with SO SOFT* inner lining
- Setting a new benchmark with POWERGUARD* protection



25867

Code	Description	Packaging
25867	FLUIDSHIELD* Level 1 Procedure Mask with SO SOFT* Lining and Earloops, Yellow	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500



25868

25868	FLUIDSHIELD* Level 1 Procedure Mask, Pleat-Style with SO SOFT* Ear Loops, Lavender	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500
-------	--	---



28802

28802	FLUIDSHIELD* Level 1 Surgical Mask with SO SOFT* Lining	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500
-------	---	---



41802

41802	FLUIDSHIELD* Level 1 Procedure Mask with Fog-Free Foam Strip and SO SOFT* Lining and Ear Loops, White	Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500
-------	---	---

MASK DISPENSERS



HALYARD* MASK DISPENSERS

Code	Description	Packaging
4966	Wall Bracket, 170mm (L) x 100mm (W) x 160mm (H)	Eaches = 1
4967	Wall Bracket, 170mm (L) x 125mm (W) x 160mm (H)	Eaches = 1
4968	Wall Bracket, 170mm (L) x 190mm (W) x 160mm (H)	Eaches = 1

TEST METHODS FOR FACE MASKS

BACTERIAL FILTRATION EFFICIENCY (BFE) AT 3 μ m: ASTM F2101

Determines the bacterial filtration efficiency (BFE) of various filtration materials. This test employs a ratio of the bacterial challenge counts to sample effluent counts to determine the percent bacterial filtration efficiency (%BFE). The test procedure allows a reproducible bacterial challenge to be delivered to test materials, and provides a standard procedure for comparison of filtration materials.

The higher the BFE, the better the filtration efficiency, the better the protection.

DIFFERENTIAL PRESSURE OR BREATHABILITY (DELTA P): 4.4.1.2 of MIL-M-36954C

Determines the air exchange differential of porous materials measured in mm H₂O/cm². This technique employs a water manometer differential upstream and downstream of the test material at a constant flow rate.

The lower the delta P, the more breathable the material is. However, this can also mean lower levels of protection.

FLUID RESISTANCE: ASTM F 1862

Standard Test Method for Resistance of Medical Face Masks to Penetration by Synthetic Blood (Horizontal Projection of Fixed Volume at Known Velocity).

Developed to simulate a challenge to the fluid splash resistance of a face mask under conditions similar to actual use. During this test, face masks are insulated with synthetic blood, and graded as a pass or fail. The masks are evaluated at three pressures: 80, 120, and 160mmHg.

The higher the pressure at which a mask passes, the greater the fluid resistance.

NATIONAL INSTITUTE FOR OCCUPATIONAL SAFETY AND HEALTH (NIOSH) Certification for Particulate Respirators: 42 CFR 84 Certified in USA

NIOSH Certification ensures the respirator has undergone the stringent testing and quality assurance requirements to meet an expected level of protection.

If a respirator has been NIOSH-certified, it will have the following markings on the respirator's packaging, user instruction insert, and/or on the respirator itself:

- The NIOSH name or logo;
- The NIOSH filter classification – in this case "N95";
- The NIOSH test and certification – or TC – approval number;
- The Approval Holder's business name (usually the manufacturer), registered trademark, or an easily understood abbreviation;
- The model or part number; and
- The manufacturing lot number.

The respirators certified under the NIOSH 42 CFR 84 regulation are tested under much more demanding conditions, providing a higher level of worker protection.



HALYARD GUARDALL* FACE SHIELD

- Designed for procedures and environments where risk of exposure to blood and/or bodily fluids through splashing is possible
- Equipped with a foam band and headband
- Can be worn with or without glasses or goggles
- Resistant to fogging

Code	Description	Packaging
41204	GUARDALL SHIELD* Full-Length Visor Only, (No Mask), Clear	Eaches/Case = 40
41205	GUARDALL SHIELD* Three-Quarter Length Visor Only (No Mask), Clear	Eaches/Case = 40



HALYARD SAFEVIEW* GLASSES

- For use in situations which may involve splashing or spraying
- Designed to be lightweight for added comfort in use
- Shipped with colorful frames to appeal to a variety of preferences
- Single-use lenses may be easily replaced into reusable frames
- Designed to be easily worn over most prescription glasses
- Constructed from natural rubber latex-free materials



Code	Description	Packaging
SV50A	SAFEVIEW* Assembled Glasses	5 Pkg. of 10 Each Case = 50 Assembled Glasses
SV100F	SAFEVIEW* Frames / Assorted Colours	10 Pkg. of 10 Each Case = 100 Frames
SV250L	SAFEVIEW* Replacement Lenses	10 Pkg. of 25 Each Case = 250 Lenses

THE WORRY-FREE CHOICE FOR EVERYONE, EVERYWHERE

Your hands are essential to the care you provide to your patients. With 5–15% of healthcare workers susceptible to latex allergy, caring for you and your hands by reducing the exposure to latex reduces the risks and harmful effects of an allergic reaction.²⁴

Halyard has created a patented, natural rubber latex-free glove to provide you with the comfort, fit and sensitivity of latex, without the added risk of allergy exposure.

It's the one glove for all your non-surgical needs, facility-wide.²⁵



INNOVATION IN PACKAGING TO REDUCE CROSS CONTAMINATION AND WASTE

As part of our commitment to reducing the risk of HAIs, Halyard introduced the SmartPULL* glove dispensing system to our STERLING* and LAVENDER* Nitrile Exam Glove boxes.

Our advanced packing technology has also reduced the waste in storage space with the larger, 200–250 count per box providing you with more for less.



The smaller tab reduces the risk of contamination and waste by limiting the contact and potential fall-out from adjacent gloves.



The larger tab enables easier access as the box empties.

NITRILE EXAM



HALYARD STERLING* NITRILE EXAM GLOVES

- For Low to High Fluid Exposure
- Ideal for Facility-Wide Use
- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (242mm)
- Enhanced Wet and Dry Grip
- Tested Against 52 Chemotherapy Agents



Code	Description	Packaging
13938	STERLING* Nitrile Exam Gloves, X-Small	Eaches/Box = 200 Boxes/Case = 10
13940	STERLING* Nitrile Exam Gloves, Small	Eaches/Box = 200 Boxes/Case = 10
13941	STERLING* Nitrile Exam Gloves, Medium	Eaches/Box = 200 Boxes/Case = 10
13942	STERLING* Nitrile Exam Gloves, Large	Eaches/Box = 200 Boxes/Case = 10
13943	STERLING* Nitrile Exam Gloves, X-Large	Eaches/Box = 170 Boxes/Case = 10



HALYARD STERLING* NITRILE-XTRA EXAM GLOVES

- For Low to High Fluid Exposure
- Ideal for Facility-Wide Use
- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 12" (310mm)
- Enhanced Wet and Dry Grip
- Tested Against 52 Chemotherapy Agents



Code	Description	Packaging
53137	STERLING* NITRILE-XTRA Exam Gloves, X-Small	Eaches/Box = 100 Boxes/Case = 10
53138	STERLING* NITRILE-XTRA Exam Gloves, Small	Eaches/Box = 100 Boxes/Case = 10
53139	STERLING* NITRILE-XTRA Exam Gloves, Medium	Eaches/Box = 100 Boxes/Case = 10
53140	STERLING* NITRILE-XTRA Exam Gloves, Large	Eaches/Box = 100 Boxes/Case = 10
53141	STERLING* NITRILE-XTRA Exam Gloves, X-Large	Eaches/Box = 100 Boxes/Case = 10

NITRILE EXAM

HALYARD PURPLE NITRILE*
EXAM GLOVES

- For Moderate to High Fluid Exposure
- Ideal for Oncology, Pharmacy, ED, Burns Units
- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (242mm)
- Enhanced Wet and Dry Grip
- Tested Against 29 Chemotherapy and other Hazardous Agents



Code	Description	Packaging
52000M	PURPLE NITRILE* Exam Gloves, X-Small	Eaches/Box = 100 Boxes/Case = 10
52001M	PURPLE NITRILE* Exam Gloves, Small	Eaches/Box = 100 Boxes/Case = 10
52002M	PURPLE NITRILE* Exam Gloves, Medium	Eaches/Box = 100 Boxes/Case = 10
52003M	PURPLE NITRILE* Exam Gloves, Large	Eaches/Box = 100 Boxes/Case = 10
52004M	PURPLE NITRILE* Exam Gloves, X-Large	Eaches/Box = 90 Boxes/Case = 10

HALYARD PURPLE NITRILE-XTRA*
HIGH RISK EXAM GLOVES

- For Moderate to High Fluid Exposure
- Ideal for Oncology, Pharmacy, ED, Burns Units
- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 12" (310mm)
- Longer Cuff for Extra Protection
- Enhanced Wet and Dry Grip
- Tested Against 52 Chemotherapy and other Hazardous Agents



Code	Description	Packaging
50600M	PURPLE NITRILE-XTRA* Exam Gloves, X-Small	Eaches/Box = 50 Boxes/Case = 10
50601M	PURPLE NITRILE-XTRA* Exam Gloves, Small	Eaches/Box = 50 Boxes/Case = 10
50602M	PURPLE NITRILE-XTRA* Exam Gloves, Medium	Eaches/Box = 50 Boxes/Case = 10
50603M	PURPLE NITRILE-XTRA* Exam Gloves, Large	Eaches/Box = 50 Boxes/Case = 10
50604M	PURPLE NITRILE-XTRA* Exam Gloves, X-Large	Eaches/Box = 50 Boxes/Case = 10

NITRILE EXAM



HALYARD PURPLE NITRILE* MAX EXAM GLOVES

- For High Fluid Exposure
- Ideal for CSSD
- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Longer length: 16" (415mm)
- Enhanced Wet and Dry Grip
- Tested Against 52 Chemotherapy and other Hazardous Agents



Code	Description	Packaging
44992	PURPLE NITRILE* MAX Exam Gloves, Small	Eaches/Box = 50 Boxes/Case = 8
44993	PURPLE NITRILE* MAX Exam Gloves, Medium	Eaches/Box = 50 Boxes/Case = 8
44994	PURPLE NITRILE* MAX Exam Gloves, Large	Eaches/Box = 50 Boxes/Case = 8
44995	PURPLE NITRILE* MAX Exam Gloves, X-Large	Eaches/Box = 50 Boxes/Case = 8



HALYARD LAVENDER* NITRILE EXAM GLOVES

- For Low to Moderate Fluid Exposure
- Ideal for Oral Health and General Use
- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (242mm)
- Enhanced Wet and Dry Grip
- Tested Against 43 Chemotherapy Agents

Code	Description	Packaging
52816	LAVENDER* Nitrile Exam Gloves, X-Small	Eaches/Box = 250 Boxes/Case = 10
52817	LAVENDER* Nitrile Exam Gloves, Small	Eaches/Box = 250 Boxes/Case = 10
52818	LAVENDER* Nitrile Exam Gloves, Medium	Eaches/Box = 250 Boxes/Case = 10
52819	LAVENDER* Nitrile Exam Gloves, Large	Eaches/Box = 250 Boxes/Case = 10
52820	LAVENDER* Nitrile Exam Gloves, X-Large	Eaches/Box = 230 Boxes/Case = 10

GLOVE DISPENSERS



50201



50202



50203

EVERYDAY GLOVE DISPENSERS

Code	Description	Packaging
50201	Exam Glove Dispenser, Clear Plastic, Single, 96mm (D) x 144mm (W) x 253mm (H)	Eaches = 1
50202	Exam Glove Dispenser, Clear Plastic, Double, 96mm (D) x 274mm (W) x 256mm (H)	Eaches = 1
50203	Exam Glove Dispenser, Clear Plastic, Triple, 96mm (D) x 406mm (W) x 256mm (H)	Eaches = 1



47467



47466

HALYARD* PURPLE NITRILE* MAX GLOVE DISPENSERS

Code	Description	Packaging
47467	Exam Glove Dispenser Extra Long (Suits 16" gloves),, Clear Plastic, Single, 75mm (D) x 139mm (W) x 409mm (H)	Eaches = 1
47466	Exam Glove Dispenser Extra Long (Suits 16" gloves),, Clear Plastic, Triple, 69 mm (D) x 399mm (W) x 409 mm (H)	Eaches = 1

BOUFFANT CAPS

- No potential fluid contact
- 1-Layer spunbond fabric
- Regular size



Code	Description	Packaging
3501	Bouffant Cap, Spunbond Fabric, Blue	Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500
3502	Bouffant Cap, Spunbond Fabric, White	Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500

SURGICAL CAP

- For light fluid contact
- Tri-Layer SMS fabric on sides
- Spunbond fabric on top
- Universal size



Code	Description	Packaging
69520	Surgical Cap, Blue	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300

SURGICAL HOOD

- For light fluid contact
- Tri-layer SMS fabric on sides
- Spunbond fabric on top
- Tie neck
- Universal size



Code	Description	Packaging
69110	Protective Surgical Hood, Blue	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300



HALYARD* OVERSHOES

Code	Description	Packaging
69252	X-TRA TRACTION* Shoe Cover, Universal	Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300
69254	X-TRA TRACTION* Shoe Cover, X-Large E	Eaches/Dispenser = 80 Dispensers/Case = 3 Eaches/Case = 240



HALYARD ANKLE GUARD* SHOE COVER

- Moderate fluid contamination
- Over the ankle coverage
- Protective tri-layer SMS fabric
- Adhesive traction strips
- Natural rubber latex-free

Code	Description	Packaging
69353	ANKLE GUARD* Shoe Cover, Tri-Layer Fabric, Adhesive Traction Strips, X-Large	Eaches/Dispenser = 90 Dispensers/Case = 3 Eaches/Case = 270




HALYARD HI-GUARD* ULTRA FULL-COVERAGE BOOT

- Potential for heavy fluid contamination
- Protection for staff when assisting with showering
- Knee high boot
- Plastic film coating over tri-layer fabric
- Two 1", 2.5cm adhesive foam traction strips
- Natural rubber latex-free

Code	Description	Packaging
69572	HI-GUARD* Full Coverage Boot, Medium	Eaches/ Dispenser = 30 Dispensers/Case = 4 Eaches/Case = 120
69672	HI-GUARD* Full Coverage Boot, Large	Eaches/ Dispenser = 30 Dispensers/Case = 4 Eaches/Case = 120



DAY SURGERY PRODUCTS
SURFACE, STAFF & PATIENT CARE





HALYARD VERSA* TOWEL

- Soft white 2 ply crepe wadding towel provides excellent absorbency and strength²⁶
- Can be used across a wide variety of medical and non-medical applications
- Wall or Bed Mounted Dispensers available for easy, economical dispensing
- ARTG 332698

Code	Size	Description	Packaging
4210	24.5cm x 41.5cm	VERSA* Towel Small	Rolls/Case = 16
4220	49cm x 41.5cm	VERSA* Towel Large	Rolls/Case = 8
7041	27.8cm W x 14.8cm H x 14.5cm D	White Plastic Dispenser Suits Code 4210	1
7056	52.5cm W x 14.8cm H x 14.5cm D	White Plastic Dispenser Suits Code 4220	1



HALYARD TERI* WIPERS

- Soft 4 ply wiper with tough polyester scrim that provide excellent absorbency and strength even when wet³⁰
- Thick 4 ply tissue provides a substantial barrier between your hands and contaminants
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698

Code	Size	Description	Packaging
4465	31.5cm x 34cm	TERI* Wipers Small	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600
4467	31.5cm x 34cm	TERI* Wipers Small Bulk Pack	Eaches/Pack = 450
4466	31.5cm x 61cm	TERI* Wipers Large	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600



HALYARD TERI* HYDRO ABSORBENT WIPERS

- Highly absorbent material – absorbs up to six times its own weight in water²⁷
- Thick cloth-like feel provides a substantial barrier between your hands and contaminants²⁹
- Ultra-low-linting³⁰
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698

Code	Size	Description	Packaging
4476	32.5cm x 33cm	TERI* Hydro Absorbent Wipers Small	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4478	32.5cm x 60cm	TERI* Hydro Absorbent Wipers Large	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450

HALYARD ISOWIPE* BACTERICIDAL WIPES



- Contains 70% Isopropyl Alcohol to provide hospital grade disinfection of clean surfaces
- ISOWIPE* Bactericidal Wipe is entered in the ARTG, 342811
- Kills COVID-19 and SARS virus in 60 seconds
- Kills Bacteria (S. aureus, E. coli, Ps. Aeruginosa, Pr. vulgaris)
- Kills Fungi (Yeast and Mould) in 10 minutes

Always read the label and follow the directions for use.



Code	Size	Description	Packaging
6835	42cm x 14cm	ISOWIPE* Bactericidal Wipes Canister	Eaches/Canister = 75 Canisters/Case = 12 Eaches/Case = 900
6836	42cm x 14cm	ISOWIPE* Refill Pack	Eaches/Refill = 75 Refills/Case = 12 Eaches/Case = 900
6837	21cm x 14cm	ISOWIPE* Minis Canister	Eaches/Canister = 150 Canisters/Case = 12 Eaches/Case = 1800
6838	21cm x 14cm	ISOWIPE* Minis Refill Pack	Eaches/Refill = 150 Refills/Case = 12 Eaches/Case = 1800
4989	17.8cm H x 9.4cm D	Stainless Steel Wall Bracket	1



HALYARD* UNDERPADS

- Soft top layer with highly absorbent layers of tissue and a strong polyethylene backing that is quiet, drapable and helps prevent strike-through of liquid
- Available in two absorbency levels – 5 ply and 8 ply – for light to moderate fluid situations
- Ideal for chair and mattress protection and as an underlay for medical procedures
- Hygienic protection that can assist in reducing the time taken for staff to clean and change linen

Code	Size	Description	Packaging
2765	39.5cm x 55cm	HALYARD* Underpads 5 Ply, Folded	Eaches/Case = 300
2768	39.5cm x 55cm	HALYARD* Underpads 8 Ply, Folded	Eaches/Case = 200



HALYARD* BED PROTECTOR

- 4.9L absorbent capacity²⁷ - ideal for moderate to heavy fluid situations
- Unique four layer design with highly absorbent core containing Super Absorbent Material, that draws fluid away from the body and locks wetness away to help keep patients dry²⁸
- Strong double layer poly base that helps prevent strikethrough²⁹ and maintain structural integrity when saturated, to keep surfaces dry and hygienic

Code	Size	Description	Packaging
4261	81cm x 78cm	HALYARD* Bed Protector	Eaches/Case = 25



HALYARD* BEDSHEET

- Soft yet thick and durable
- Strong and hygienic cover
- Continuous roll
- Wall or bed mounted dispenser available for economical dispensing

Code	Size	Description	Packaging
4260	54cm x 80m	HALYARD* Bedsheet Roll	Rolls/Case = 6
4924		Dispenser White Steel Suits Code 4260	1



HALYARD* BENCHROLL

- Ideal for protecting bench and surfaces from staining and other spills
- Highly absorbent layers of tissue and one layer of strong polyethylene backing to help prevent strike-through of liquids
- Embossed pattern increases liquid absorbency so spills are quickly soaked up, helping to prevent run off
- Economical – cut to desired length to reduce unnecessary waste

Code	Size	Description	Packaging
2911	41.5cm x 91m	HALYARD* Benchroll 4 Ply Roll	Rolls/Case = 2



HALYARD* SINGLE USE LINEN

Code	Size	Description	Packaging
67783	101.5cm x 213cm (40" x 84"), Bulk	Stretcher Sheet, Blue	80 eaches/case
67843	101.5cm x 183cm (40" x 72"), Bulk	Draw Sheet, Poly-Coated* Blue *Coated with Polyethylene Film	80 eaches/case
67803	50.5cm x 73.5cm (20in x 29in), Bulk	Pillowcase, Blue	200 eaches/case



HALYARD* NEUTRAL HAND WIPES

- Ideal for limited mobility patients
- Alcohol-free wipes containing neutral ingredients to clean visibly soiled hands when access to soap and water is limited
- Contains Vitamin E to restore moisture to skin for frequently washed hands
- Resealable pack
- ARTG 273333

Code	Size	Description	Packaging
8835	20cm x 20cm	HALYARD* Neutral Hand Wipes	Eaches/Pack = 15 Packs/Case = 40 Eaches/Case = 600



HALYARD VERSA* TOWEL

- White 2 ply crepe towel provides excellent absorbency and strength²²
- Soft bulky crepe towel is gentle on frequently washed hands
- Wall or Bed Mounted Dispensers available for easy, economical dispensing
- ARTG 332698

Code	Size	Description	Packaging
4210	24.5cm x 41.5cm	VERSA* Towel Small	Rolls/Case = 16
4220	49cm x 41.5cm	VERSA* Towel Large	Rolls/Case = 8
7041	27.8cm W x 14.8cm H x 14.5cm D	White Plastic Dispenser Suits Code 4210	1
7056	52.5cm W x 14.8cm H x 14.5cm D	White Plastic Dispenser Suits Code 4220	1



HALYARD TERI* WIPERS

- Thick, soft 4 ply tissue – gentle on frequently washed hands
- Excellent wet strength with scrim reinforcement²²
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698

Code	Size	Description	Packaging
4465	31.5cm x 34cm	TERI* Wipers Small	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600
4467	31.5cm x 34cm	TERI* Wipers Small Bulk Pack	Eaches/Pack = 450
4466	31.5cm x 61cm	TERI* Wipers Large	Eaches/Pack = 100 Packs/Case = 6 Eaches/Case = 600



HALYARD TERI* HYDRO ABSORBENT WIPERS

- Highly absorbent material – absorbs up to six times its own weight in water²³
- Ultra low-linting²⁶
- Thick cloth-like feel is gentle on frequently washed hands
- Perforated box for easy dispensing and the convenience of portability
- ARTG 332698

Code	Size	Description	Packaging
4476	32.5cm x 33cm	TERI* Hydro Absorbent Wipers Small	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450
4478	32.5cm x 60cm	TERI* Hydro Absorbent Wipers Large	Eaches/Pack = 75 Packs/Case = 6 Eaches/Case = 450

NOTE: Please see package inserts for directions for use for each of these products.
Always read the label and follow the directions for use.



HALYARD* WASHCLOTH EXTRA

- Ideal as patient face or body washer
- Soft, thick and super absorbent - absorbs up to eight times its weight in water²³
- Strong and durable enough to not tear easily²² in wet or dry applications
- Low linting²⁶
- Clear resealable bag to help prevent dust and other particles from contaminating cloths

Code	Size	Description	Packaging
2711	35cm x 40cm	HALYARD* Washcloth Extra	Eaches/Pack = 20 Packs/Case = 30 Eaches/Case = 600



HALYARD* PROTECTAPAD

- Ideal for use as a bib to protect patients' clothing and hygienic protection for patient surfaces
- Highly absorbent layers of tissue and one layer of strong poly backing to help prevent strike-through of liquids
- Embossed pattern increases liquid absorbency so spills are quickly soaked up, helping to prevent run off
- Soft and drapable for optimal patient comfort and use on varied surfaces

Code	Size	Description	Packaging
2706	28.5cm x 21.5cm	HALYARD* Protectapad Small	Eaches/Pack = 100 Packs/Case = 8 Eaches/Case = 800
2705	28.5cm x 43cm	HALYARD* Protectapad Large	Eaches/Pack = 100 Packs/Case = 4 Eaches/Case = 400



HALYARD* BILATERAL FACIAL ICE PACK

- 2 ice packs bonded together to evenly deliver cold therapy to both sides of the jaw
- Can be used for recovery after maxillofacial and dental procedures
- Single-patient use

Code	Description	Size	Packaging
33101	Bilateral Facial Ice Pack	13cm x 30cm	12 eaches/box 2 boxes/case 24 eaches/case



HALYARD* EYECARE (EENT) ICE PACK

- Designed for easy application after plastic surgery, ophthalmic and sinus procedures
- Adjustable strap ensures the pack sits comfortably and securely on the face
- Single-patient use

Code	Description	Size	Packaging
33150	EyeCare (EENT) Ice Pack, 1 Strap	11.4cm x 25cm	15 eaches/box 2 boxes/case 30 eaches/case



HALYARD* STAY-DRY ICE PACK

- For general use
- Delivers up to 2.5 hours of consistent, dry cold therapy solutions
- STAY-DRY construction features 3 layers which help to reduce wetness due to condensation
- 4 ties for versatile application on shoulders, knees, elbows, etc
- Single-patient use

Code	Description	Size	Packaging
33495	STAY-DRY Ice Pack, Small, 4 Ties	13cm x 30cm	25 eaches/box 2 boxes/case 50 eaches/case
33500	STAY-DRY Ice Pack, Large, 4 Ties	16.5cm x 36cm	25 eaches/box 2 boxes/case 50 eaches/case

I.V. ARTERIAL AND WRIST SUPPORTS



HALYARD* NEONATAL & INFANT/CHILD I.V. SUPPORT

- Padded for comfort
- Aluminium core easily bends to conform to various angles of small wrists and feet
- Single-patient use

Code	Description	Size	Packaging
29965	Neonatal I.V. Support	9cm x 4cm, 13cm straps	30 eaches/case
29967	Infant/Child I.V. Support	16cm x 5cm, 23cm straps	30 eaches/case



HALYARD* HAND-AID* ARTERIAL WRIST SUPPORT

- Uniquely anatomically designed to position the wrist for accurate and efficient application of arterial lines
- The instant support of the foam block saves setup time
- Single-patient use

Code	Description	Size	Packaging
29990	Hand-Aid* Arterial Wrist Support, Paediatric	16cm x 7cm, 23cm straps	20 eaches/case
29980	Hand-Aid* Arterial Wrist Support, Adult	23cm x 9cm, 37cm straps	20 eaches/case
29985	Hand-Aid* Arterial Wrist Support, Replacement Straps & Starter Block, Compatible with 29980		50 eaches/case



HALYARD* HAND-AID* WRIST SUPPORT

- Anatomically designed to position the wrist for accurate and efficient cannulisation in dorsal hand veins
- Single-patient use

Code	Description	Size	Packaging
29970	Hand-Aid* I.V. Wrist Support, Paediatric	14cm x 5cm, 23cm straps	20 eaches/case
29900	Hand-Aid* I.V. Wrist Support, Adult	16cm x 7cm, 23cm straps	20 eaches/case



HALYARD* ANTECUBITAL I.V. SUPPORT

- The unique and comfortable way to immobilise and support during antecubital I.V. therapy
- Single-patient use

Code	Description	Size	Packaging
29920	Antecubital I.V. Support	32cm x 11cm, 48 cm straps	10 eaches/case

1. Data on File.
2. Comparing like fabric weights. Data on File.
3. Physical Properties Testing Dossier H01104.
4. AORN 2012 Perioperative Standards, Recommended Practices and Guidelines.
5. Data on File.
6. Data on File.
7. Kurt Salmon & Associates, Inc. Study to quantify the time taken to wrap and open packages using Kimguard One-Step* compared to conventional sequential wrap.
8. Data on File.
9. Whole of Packaging Microbial Aerosol Challenge Test. Data on File.
10. Maintenance of Package Integrity Study: Dunkleberg, Leximed Reports, 08-L046 to 09-L024.
11. EQ-STM-4566.
12. Comparison of Smart-Fold* Wrap to One-Step* Wrap and Sequential Wrap: Kurt Salmon and Associates, Time and Motion Studies on the Labour Requirements of the Smart-Fold, One-Step and Sequential Sterilisation Wraps, June 2011.
13. Harter B. Infection Control Today 2003; May Edition
14. Dunkelberg H. and Rohmanns, S. Infection Control and Hospital Epidemiology 2006; 27(4):367-371.
15. TNO report, May 1999, Determination of the microbial barrier properties of twelve types of packaging material according to the Final Pack Method.
16. TNO report, V6412 April 2005, Determination of the microbial barrier properties of thirty types of packaging material from Kimberly-Clark, according to the Final Pack Method. TNO (Netherlands Organization for Applied Scientific Research) is a not-for-profit organisation in the Netherlands that focuses on applied science. As a statutory organisation, TNO has an independent position that allows to give objective, scientifically founded judgements.
17. AORN Journal 2000; 72(6): 1052-1059.
18. Data on file. Claim-07729.
19. Tyvek 1073B Medical Device Packaging
<https://www.dupont.com/products/tyvek-1073b-medical-device-packaging.html>
20. Billerudkorsnäs Technical Data Sheet.
21. ISO 11607, EN868, ANSI/AAMI ST79, AS/NZS 4187.
22. Per Hydrohead test method. Data on File.
23. As measured by Hydrohead. ASTM1760 Synthetic Blood Penetration. Data on File.
24. 2010. The Australasian Society of Clinical Immunology and Allergy (ASCI) <http://www.allergy.org.au/health-professionals/hp-information/asthma-and-allergy/latex-allergy> (Viewed 26/02/15)
25. Halyard offers a line of gloves specifically for use with chemotherapeutic agents.
26. Mechanical break resistance test data on file.[†]
27. Absorption capacity test data on file.[†]
28. Re-wet test data on file.[†]
29. Halyard Product Specifications data on file.
30. Gelbo lint test data on file.[†]

[†] These tests are directional based on small sample sizes.



Infection Prevention Product Solutions
Halyard* Education Foundation
Tools & Best Practice
Commitment to Excellence and Sustainability
Experienced Customer Support
Knowledgeable Sales Force
Innovation & Technology



O&M Halyard Australia Pty. Limited

Customer Care:

Phone: 1800 664 227

Email: CustomerService_ANZ@hyh.com

O&M Halyard makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

www.halyardhealth.com.au



* Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2022 O&M Halyard, Inc. All Rights Reserved.

DuPont™, the DuPont Oval Logo and Tyvek® are trademarks or registered trademarks of DuPont or its affiliates. Copyright © 2022 DuPont de Nemours Inc.

^ Registered Trademark or Trademark of Kimberly-Clark Worldwide, Inc. and under license to O&M Halyard, Inc. 22:06:WFHY947v16