



PPE CATALOGUE

AUSTRALIA &
NEW ZEALAND



QUICK LINKS

FACIAL PROTECTION

| | |
|------------------------------------|----------|
| Levels of Facial Protection | page 4 |
| Masks: Level 3 Protection: | |
| FLUIDSHIELD* N95 Face Masks | page 5 |
| FLUIDSHIELD* Face Masks | page 5 |
| Level 2 Protection: | |
| FLUIDSHIELD* Face Masks | page 6 |
| Level 1 Protection: | |
| FLUIDSHIELD* Face Masks | page 6 |
| Standard Face Masks | page 7–8 |
| Eyewear: | |
| GUARDALL* Face Shields | page 9 |
| SAFEVIEW* Glasses | page 9 |
| Accessories | page 10 |
| Test Methods for Face Masks | page 11 |

EXAMINATION GLOVES

| | |
|------------------------------|---------|
| The Worry-Free Choice | page 13 |
| STERLING* Nitrile | page 14 |
| STERLING* Nitrile-Xtra | page 14 |
| PURPLE NITRILE* | page 15 |
| PURPLE NITRILE-XTRA* | page 15 |
| PURPLE NITRILE* MAX | page 16 |
| LAVENDER* Nitrile | page 16 |
| AQUASOFT* Nitrile | page 17 |
| Glove Dispensers | page 18 |

APPAREL

| | |
|---|---------|
| Sequence for Donning and Doffing Personal Protective Equipment (PPE) | page 20 |
| Gowns: | |
| Isolation Gowns | page 21 |
| Impervious Gowns | page 22 |
| Chemotherapy Gowns | page 23 |
| Scrub Suits | page 24 |
| Headwear | page 25 |
| Footwear | page 26 |



HALYARD, YOUR PARTNER IN COMPLIANCE TO THE MASK STANDARD

AS 4381:2015 IS THE STANDARD FOR SINGLE-USE FACE MASKS

Studies show that fluid strikes the face area of OR staff an average 45–51% of the time.¹ Surgeons are particularly at risk, with 90% not aware they've been splashed.²

This puts you, your staff and your patients at considerable risk.


NEW STANDARD AS 4381:2015 (RED INDICATES CHANGE)

| CHARACTERISTICS | LEVEL 1 1 | LEVEL 2 2 | LEVEL 3 3 | TEST METHOD |
|---|--|---|--|-------------------------------------|
| | Level 1 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the minimum velocity specified in Table 2, bacterial filtration efficiency and differential pressure. APPLICATIONS: For general purpose medical procedures, where the wearer is not at risk of blood or bodily fluid splash or to protect staff and/or the patient from droplet exposure to microorganisms (e.g. patient with upper respiratory tract infection visits GP) | Level 2 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the middle velocity specified in Table 2, bacterial filtration efficiency and differential pressure. APPLICATIONS: For use in emergency departments, dentistry, changing dressings on small or healing wounds where minimal blood droplet exposure may possibly occur (e.g. endoscopy procedures) | Level 3 barrier medical face mask materials are evaluated for resistance to penetration by synthetic blood at the maximum velocity specified in Table 2, bacterial filtration efficiency and differential pressure. APPLICATIONS: For all surgical procedures, major trauma first aid or in any area where the health care worker is at risk of blood or bodily fluid splash (e.g. orthopaedic, cardiovascular procedures) | |
| Bacterial Filtration Efficiency (BFE) % | ≥ 95% | ≥ 98% | ≥ 98% | ASTM F2101-14 or EN 14683:2014 |
| Particulate Filtration Efficiency (PFE) % (0.1 µm) | Not Required | Not Required | Not Required | N/A |
| Differential Pressure (Delta P) mm H ₂ O/cm ² | < 4.0 | < 5.0 | < 5.0 | EN 14683:2014 |
| Resistance to penetration by synthetic blood (fluid resistance) min pressure in mm Hg for pass result | 80mm Hg | 120mm Hg | 160mm Hg | ASTM F1862 / F1862M-13 or ISO 22669 |

YOU, YOUR STAFF AND YOUR PATIENTS NEED THE RIGHT FACIAL PROTECTION


Halyard face masks give you the confidence that you and your patients are protected.

- ✓ All HALYARD* FLUIDSHIELD* Masks meet the new Australian Standard AS 4381:2015.
- ✓ FLUIDSHIELD* Masks are available in Level 1, Level 2 and Level 3.
- ✓ The level of protection is clearly labelled on the mask nose piece and on the box, so you and your staff can easily identify the right mask for each task.



SUPERIOR PROTECTION NEVER FELT SO GOOD

FLUIDSHIELD* masks not only provide superior cover, with four layers of barrier protection; they also feature an inner SO SOFT* lining to deliver the critical balance between protection and clinician comfort.




1. & 2. CG Davies et al, Ann R Coll Surg Engl 2007 89(8): 770–772; S Endo et al, Jour Hospital Infection 2007 (67) 1: 56–61; DC Berridge et al, Br J Surg 1993 80 (11): 1379–80.

FACIAL PROTECTION

For the highest levels of protection, make sure you choose HALYARD* Face Masks and Eyewear.

NIOSH APPROVED N95 RESPIRATOR FACE MASKS



- Provides fluid resistance at 160mm Hg
- Tested against 7 chemotherapy drugs
- The highest level of protection available
- Recommended for Laser Procedures



46727

| Code | Description | Packaging |
|-------|--|--|
| 46727 | FLUIDSHIELD* Level 3 N95 Particulate Filter Respirator and Surgical Mask, Pouch-Style with Headbands, Regular Size, Orange | Eaches/Box = 35 Boxes/Case = 6 Eaches/Case = 210 |
| 46827 | FLUIDSHIELD* Level 3 N95 Particulate Filter Respirator and Surgical Mask, Pouch-Style with Headbands, Small Size, Orange | Eaches/Box = 35 Boxes/Case = 6 Eaches/Case = 210 |



FLUIDSHIELD* FACE MASKS

- Provides fluid resistance at 160mm Hg
- Level 3 printed on mask nose piece and box
- Four layers of fluid resistance and bacterial challenge
- Superior comfort without compromising your protection, with SO SOFT* inner lining
- Setting a new benchmark with POWERGUARD* protection



28797



28800



28804



47107



47147



48207



48207



48297

| Code | Description | Packaging |
|-------|---|---|
| 28797 | FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with SO SOFT* Ear Loops | Eaches/Box = 40 Boxes/Case = 10 Eaches/Case = 400 |
| 28800 | FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Wraparound Visor and SO SOFT* Lining and SO SOFT* Ear Loops | Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100 |
| 28804 | FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with Wraparound Visor and SO SOFT* Lining | Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100 |
| 47107 | FLUIDSHIELD* Level 3 Fog-Free Procedure Mask, Pleat-Style with SO SOFT* Ear Loops, Foam Band, Orange | Eaches/Box = 40 Boxes/Case = 10 Eaches/Case = 400 |
| 47147 | FLUIDSHIELD* Level 3 Fog-Free Procedure Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with SO SOFT* Ear Loops, Foam Band, Orange | Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100 |
| 48207 | FLUIDSHIELD* Level 3 Fog-Free Surgical Mask, Pleat-Style with Ties, Orange | Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300 |
| 48247 | FLUIDSHIELD* Level 3 Fog-Free Surgical Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Ties, Orange | Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100 |
| 48297 | FLUIDSHIELD* Level 3 Care Bear Surgical Mask, Pleat-Style with Ties, Pink and Blue Teddy Bear Print | Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300 |

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

FACIAL PROTECTION



FLUIDSHIELD* FACE MASKS

- Provides fluid resistance at 120mm Hg
- Level 2 printed on mask nose piece and box
- Four layers of fluid resistance and bacterial challenge
- Superior comfort without compromising your protection, with SO SOFT* inner lining
- Setting a new benchmark with POWERGUARD* protection



62114

| Code | Description | Packaging |
|-------|---|--|
| 62114 | FLUIDSHIELD* Level 2 Fog-Free Surgical Mask with Wraparound SPLASHGUARD* Visor, Pleat-Style with Ties, Blue with Orange Diamond Print | Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100 |



62115

| | | |
|-------|--|--|
| 62115 | FLUIDSHIELD* Level 2 Fog-Free Procedure Mask with SO SOFT* Ear Loops, Blue | Eaches/Box = 25 Boxes/Case = 4 Eaches/Case = 100 |
|-------|--|--|



FLUIDSHIELD* FACE MASKS

- Provides fluid resistance at 80mm Hg
- Level 1 printed on mask nose piece and box
- Four layers of fluid resistance and bacterial challenge
- Superior comfort without compromising your protection with SO SOFT* inner lining
- Setting a new benchmark with POWERGUARD* protection



25868

| Code | Description | Packaging |
|-------|--|---|
| 25868 | FLUIDSHIELD* Level 1 Procedure Mask, Pleat-Style with SO SOFT* Ear Loops, Lavender | Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500 |



28802

| | | |
|-------|---|---|
| 28802 | FLUIDSHIELD* Level 1 Surgical Mask with SO SOFT* Lining | Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500 |
|-------|---|---|



25867

| | | |
|-------|--|---|
| 25867 | FLUIDSHIELD* Level 1 Procedure Mask with SO SOFT* Lining & Ear Loops, Yellow | Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500 |
|-------|--|---|



25867

| | | |
|-------|--|---|
| 41802 | FLUIDSHIELD* Level 1 Procedure Mask with Fog-Free Foam Strip & SO SOFT* Lining, Ear Loops, White | Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500 |
|-------|--|---|

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

FACIAL PROTECTION

STANDARD FACE MASKS

HALYARD* SURGICAL MASKS

- Protective three-layer mask construction, for when exposure to blood and/or bodily fluids is not a risk
- Appropriate for general patient care and isolation
- Tie attachment



47265

| Code | Description | Packaging |
|-------|---|--|
| 47625 | LAZER* Surgical Mask, Pleat-Style With Ties, Silver | Eaches/Box = 35 Boxes/Case = 6 Eaches/Case = 210 |



48100

| | | |
|-------|--|--|
| 48100 | THE LITE ONE* Surgical Mask, Pleat-Style with Ties, Blue | Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300 |
|-------|--|--|



37525

| | | |
|-------|--|--|
| 37525 | Surgical Mask, Expanded Chamber, Green | Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300 |
|-------|--|--|



48390

| | | |
|-------|--|--|
| 48390 | SO SOFT* Surgical Mask, Pleat-Style with Ties, White | Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300 |
|-------|--|--|



49214

| | | |
|-------|---|--|
| 49214 | Fog-Free Surgical Mask, Pleat-Style with Ties, Blue | Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300 |
|-------|---|--|



49215

| | | |
|-------|--|--|
| 49215 | Fog-Free Surgical Mask with DERMA-TOUCH* Tape, Film Strip, Green | Eaches/Box = 50 Boxes/Case = 6 Eaches/Case = 300 |
|-------|--|--|

FACIAL PROTECTION

HALYARD* PROCEDURE MASKS

- Protective three-layer mask construction, for when exposure to blood and/or bodily fluids is not a risk
- Appropriate for general patient care and isolation
- Earloop attachment



47080



47127



62356

| Code | Description | Packaging |
|-------|--|---|
| 47080 | Procedure Mask, Pleat-Style with Ear Loops, Blue | Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500 |
| 47127 | Disney® Child's Face Mask, Pleat-Syle with Ear Loops, Multi-Colour | Eaches/Box = 75 Boxes/Case = 10 Eaches/Case = 750 |
| 62356 | THE LITE ONE* Procedure Mask, Pleat-Style with Ear Loops, Blue | Eaches/Box = 50 Boxes/Case = 10 Eaches/Case = 500 |

FACIAL PROTECTION

EYEWEAR

GUARDALL* FACE SHIELD

- Designed for procedures and environments where risk of exposure to blood and/or bodily fluids through splashing is possible
- Equipped with a foam band and headband
- Can be worn with or without glasses or goggles
- Resistant to fogging



| Code | Description | Packaging |
|-------------|---|------------------|
| 41204 | GUARDALL SHIELD* Full-Length Visor Only, (No Mask), Clear | Eaches/Case = 40 |
| 41205 | GUARDALL SHIELD* Three-Quarter Length Visor Only (No Mask), Clear | Eaches/Case = 40 |

SAFEVIEW* GLASSES

- For use in situations which may involve splashing or spraying
- Designed to be lightweight for added comfort in use
- Shipped with colorful frames to appeal to a variety of preferences
- Single-use lenses may be easily replaced into reusable frames
- Designed to be easily worn over most prescription glasses
- Constructed from natural rubber latex-free materials



| Code | Description | Packaging |
|-------------|-------------------------------------|--|
| SV50A | SAFEVIEW* Assembled Glasses | 5 Pkg. of 10 Each Case = 50 Assembled Glasses |
| SV100F | SAFEVIEW* Frames / Assorted Colours | 10 Pkg. of 10 Each Case = 100 Frames |
| SV250L | SAFEVIEW* Replacement Lenses | 10 Pkg. of 25 Each Case = 250 Lenses |

FACIAL PROTECTION

ACCESSORIES

HALYARD* N95 RESPIRATOR ACCESSORIES

- Contain Instructions and Materials to Conduct Fit Testing for HALYARD PFR95* N95 Respirators



47950



99009



99009

| Code | Description | Packaging |
|-------|--------------------------|----------------------------------|
| 47950 | Qualitative Fit Test Kit | Eaches/Kit = 1 |
| 99003 | Saccharin Solutions | Eaches/Box = 2 Boxes/Case = 1 |
| 99009 | BITREX* Solutions | Eaches/Box = 2 Boxes/Case = 1 |
| 99018 | Nebulisers | Eaches/Box = 2 Boxes/Case = 1 |

HALYARD* MASK DISPENSER WALL BRACKET



| Code | Description | Packaging |
|------|---|------------|
| 4966 | Wall Bracket, 170mm (L) x 100mm (W) x 160mm (H) | Eaches = 1 |
| 4967 | Wall Bracket, 170mm (L) x 125mm (W) x 160mm (H) | Eaches = 1 |
| 4968 | Wall Bracket, 170mm (L) x 190mm (W) x 160mm (H) | Eaches = 1 |

FACIAL PROTECTION

TEST METHODS FOR FACE MASKS

BACTERIAL FILTRATION EFFICIENCY (BFE) AT 3µm: ASTM F2101

Determines the bacterial filtration efficiency (BFE) of various filtration materials. This test employs a ratio of the bacterial challenge counts to sample effluent counts to determine the percent bacterial filtration efficiency (%BFE). The test procedure allows a reproducible bacterial challenge to be delivered to test materials, and provides a standard procedure for comparison of filtration materials.

The higher the BFE, the better the filtration efficiency, the better the protection.

DIFFERENTIAL PRESSURE OR BREATHABILITY (DELTA P): 4.4.1.2 of MIL-M-36954C

Determines the air exchange differential of porous materials measured in mm H₂O/cm². This technique employs a water manometer differential upstream and downstream of the test material at a constant flow rate.

**The lower the delta P, the more breathable the material is.
However, this can also mean lower levels of protection.**

FLUID RESISTANCE: ASTM F 1862

Standard Test Method for Resistance of Medical Face Masks to Penetration by Synthetic Blood (Horizontal Projection of Fixed Volume at Known Velocity).

Developed to simulate a challenge to the fluid splash resistance of a face mask under conditions similar to actual use. During this test, face masks are insulted with synthetic blood, and graded as a pass or fail. The masks are evaluated at three pressures: 80, 120, and 160mmHg.

The higher the pressure at which a mask passes, the greater the fluid resistance.

NATIONAL INSTITUTE FOR OCCUPATIONAL SAFETY AND HEALTH (NIOSH) Certification for Particulate Respirators: 42 CFR 84 Certified in USA

NIOSH Certification ensures the respirator has undergone the stringent testing and quality assurance requirements to meet an expected level of protection.

If a respirator has been NIOSH-certified, it will have the following markings on the respirator's packaging, user instruction insert, and/or on the respirator itself:

- The NIOSH name or logo;
- The NIOSH filter classification – in this case "N95";
- The NIOSH test and certification – or TC – approval number;
- The Approval Holder's business name (usually the manufacturer), registered trademark, or an easily understood abbreviation;
- The model or part number; and
- The manufacturing lot number.

The respirators certified under the NIOSH 42 CFR 84 regulation are tested under much more demanding conditions, providing a higher level of worker protection.



EXAMINATION GLOVES



THE WORRY-FREE CHOICE FOR EVERYONE, EVERYWHERE

Your hands are essential to the care you provide to your patients. With 5–15% of healthcare workers susceptible to latex allergy, caring for you and your hands by reducing the exposure to latex reduces the risks and harmful effects of an allergic reaction¹.

Halyard has created a patented, natural rubber latex-free glove to provide you with the comfort, fit and sensitivity of latex, without the added risk of exposure.

It's the one glove for all your non-surgical needs, facility-wide².



INNOVATION IN PACKAGING TO REDUCE CROSS CONTAMINATION AND WASTE

As part of our commitment to reducing the risk of HAIs, Halyard introduced the SmartPULL* glove dispensing system to our STERLING* and LAVENDER* Nitrile Exam Glove boxes.

Our advanced packing technology has also reduced the waste in storage space with the larger, 200–250 count per box providing you with more for less.



The smaller tab reduces the risk of contamination and waste by limiting the contact and potential fall-out from adjacent gloves.



The larger tab enables easier access as the box empties.

REFERENCES:

- 2010, The Australasian Society of Clinical Immunology and Allergy (ASCIA) <http://www.allergy.org.au/health-professionals/hp-information/asthma-and-allergy/latex-allergy> (Viewed 26/02/15)
- Halyard offers a line of gloves specifically for use with chemotherapeutic agents.

EXAMINATION GLOVES

NITRILE EXAM GLOVES

STERLING* NITRILE EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- Silver
- CE Marked
- Enhanced Wet and Dry Grip
- Tested Against 52 Chemotherapy Agents



| Code | Description | Packaging |
|-------|--|-------------------------------------|
| 13938 | STERLING* Nitrile Exam Gloves, X-Small | Eaches/Box = 200 Boxes/Case = 10 |
| 13940 | STERLING* Nitrile Exam Gloves, Small | Eaches/Box = 200 Boxes/Case = 10 |
| 13941 | STERLING* Nitrile Exam Gloves, Medium | Eaches/Box = 200 Boxes/Case = 10 |
| 13942 | STERLING* Nitrile Exam Gloves, Large | Eaches/Box = 200 Boxes/Case = 10 |
| 13943 | STERLING* Nitrile Exam Gloves, X-Large | Eaches/Box = 170 Boxes/Case = 10 |

STERLING* NITRILE-XTRA EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 12" (305mm)
- Silver
- Enhanced Wet and Dry Grip
- Tested Against 52 Chemotherapy Agents



| Code | Description | Packaging |
|-------|---|-------------------------------------|
| 53137 | STERLING* NITRILE-XTRA Exam Gloves, X-Small | Eaches/Box = 100 Boxes/Case = 10 |
| 53138 | STERLING* NITRILE-XTRA Exam Gloves, Small | Eaches/Box = 100 Boxes/Case = 10 |
| 53139 | STERLING* NITRILE-XTRA Exam Gloves, Medium | Eaches/Box = 100 Boxes/Case = 10 |
| 53140 | STERLING* NITRILE-XTRA Exam Gloves, Large | Eaches/Box = 100 Boxes/Case = 10 |
| 53141 | STERLING* NITRILE-XTRA Exam Gloves, X-Large | Eaches/Box = 100 Boxes/Case = 10 |

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

EXAMINATION GLOVES

PURPLE NITRILE* EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- Purple
- CE Marked
- Enhanced Wet and Dry Grip
- Tested Against 29 Chemotherapy Agents



| Code | Description | Packaging |
|--------|--------------------------------------|-------------------------------------|
| 52000M | PURPLE NITRILE* Exam Gloves, X-Small | Eaches/Box = 100 Boxes/Case = 10 |
| 52001M | PURPLE NITRILE* Exam Gloves, Small | Eaches/Box = 100 Boxes/Case = 10 |
| 52002M | PURPLE NITRILE* Exam Gloves, Medium | Eaches/Box = 100 Boxes/Case = 10 |
| 52003M | PURPLE NITRILE* Exam Gloves, Large | Eaches/Box = 100 Boxes/Case = 10 |
| 52004M | PURPLE NITRILE* Exam Gloves, X-Large | Eaches/Box = 90 Boxes/Case = 10 |

PURPLE NITRILE-XTRA* HIGH RISK EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Purple
- Length: 12" (305mm)
- Longer Cuff for Extra Protection
- CE Marked
- Enhanced Wet and Dry Grip
- Tested Against 52 Chemotherapy Agents



| Code | Description | Packaging |
|--------|---|------------------------------------|
| 50600M | PURPLE NITRILE-XTRA* Exam Gloves, X-Small | Eaches/Box = 50 Boxes/Case = 10 |
| 50601M | PURPLE NITRILE-XTRA* Exam Gloves, Small | Eaches/Box = 50 Boxes/Case = 10 |
| 50602M | PURPLE NITRILE-XTRA* Exam Gloves, Medium | Eaches/Box = 50 Boxes/Case = 10 |
| 50603M | PURPLE NITRILE-XTRA* Exam Gloves, Large | Eaches/Box = 50 Boxes/Case = 10 |
| 50604M | PURPLE NITRILE-XTRA* Exam Gloves, X-Large | Eaches/Box = 50 Boxes/Case = 10 |

NOTE: Please see package inserts for directions for use for each of these products
Always read the label and follow the directions for use..

EXAMINATION GLOVES

PURPLE NITRILE* MAX EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Purple
- Longer Length: 16 inches¹ (406mm)
- CE Marked
- Enhanced Wet and Dry Grip
- Heavier weight: 20 grams²
- Increased fingertip thickness
- Tested Against 52 Chemotherapy Agents and 25 Commonly Used Chemical and Detergents



| Code | Description | Packaging |
|-------|--|-----------------------------------|
| 44992 | PURPLE NITRILE* MAX Exam Gloves, Small | Eaches/Box = 50 Boxes/Case = 8 |
| 44993 | PURPLE NITRILE* MAX Exam Gloves, Medium | Eaches/Box = 50 Boxes/Case = 8 |
| 44994 | PURPLE NITRILE* MAX Exam Gloves, Large | Eaches/Box = 50 Boxes/Case = 8 |
| 44995 | PURPLE NITRILE* MAX Exam Gloves, X-Large | Eaches/Box = 50 Boxes/Case = 8 |

LAVENDER* NITRILE EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- Lavender
- Enhanced Wet and Dry Grip
- Tested Against 43 Chemotherapy Agents



| Code | Description | Packaging |
|-------|--|-------------------------------------|
| 52816 | LAVENDER* Nitrile Exam Gloves, X-Small | Eaches/Box = 250 Boxes/Case = 10 |
| 52817 | LAVENDER* Nitrile Exam Gloves, Small | Eaches/Box = 250 Boxes/Case = 10 |
| 52818 | LAVENDER* Nitrile Exam Gloves, Medium | Eaches/Box = 250 Boxes/Case = 10 |
| 52819 | LAVENDER* Nitrile Exam Gloves, Large | Eaches/Box = 250 Boxes/Case = 10 |
| 52820 | LAVENDER* Nitrile Exam Gloves, X-Large | Eaches/Box = 230 Boxes/Case = 10 |

REFERENCES:

1. Specification target. Data on file
2. X-Large size

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

EXAMINATION GLOVES

AQUASOFT* NITRILE EXAM GLOVES

- Powder-Free
- Natural Rubber Latex-Free
- Non-Sterile
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Length: 9.5" (241mm)
- CE Marked
- Textured fingertips
- Tested Against 34 Chemotherapy Agents
- Higher count per box saves space and waste: 3000 units per case / 300 gloves per inner



| Code | Description | Packaging |
|-------|--|-------------------------------------|
| 47634 | AQUASOFT* Nitrile Exam Gloves, X-Small | Eaches/Box = 300 Boxes/Case = 10 |
| 47635 | AQUASOFT* Nitrile Exam Gloves, Small | Eaches/Box = 300 Boxes/Case = 10 |
| 47636 | AQUASOFT* Nitrile Exam Gloves, Medium | Eaches/Box = 300 Boxes/Case = 10 |
| 47637 | AQUASOFT* Nitrile Exam Gloves, Large | Eaches/Box = 300 Boxes/Case = 10 |
| 47638 | AQUASOFT* Nitrile Exam Gloves, X-Large | Eaches/Box = 250 Boxes/Case = 10 |

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.

GLOVE DISPENSERS



HALYARD* EVERYDAY GLOVE DISPENSERS

| Code | Description | Packaging |
|-------|---|------------|
| 50201 | Exam Glove Dispenser (Suits 9.5" - 12" gloves), Clear Plastic, Single, 96mm (D) x 144mm (W) x 253mm (H) | Eaches = 1 |
| 50202 | Exam Glove Dispenser (Suits 9.5" - 12" gloves), Clear Plastic, Double, 96mm (D) x 274mm (W) x 256mm (H) | Eaches = 1 |
| 50203 | Exam Glove Dispenser (Suits 9.5" - 12" gloves), Clear Plastic, Triple, 96mm (D) x 406mm (W) x 256mm (H) | Eaches = 1 |

HALYARD* PURPLE NITRILE* MAX GLOVE DISPENSERS



| Code | Description | Packaging |
|-------|---|------------|
| 47467 | Exam Glove Dispenser Extra Long (Suits 16" gloves), Clear Plastic, Single, 75mm (D) x 139mm (W) x 409mm (H) | Eaches = 1 |
| 47466 | Exam Glove Dispenser Extra Long (Suits 16" gloves), Clear Plastic, Triple, 69 mm (D) x 399mm (W) x 409 mm (H) | Eaches = 1 |



SEQUENCE FOR DONNING PERSONAL PROTECTIVE EQUIPMENT (PPE):

1. PROTECT YOUR BODY

Ensure the gown:

- Covers torso from neck to knees
- Covers arms to end of wrists
- Wraps around the sides to fasten at neck and back



2. PROTECT YOUR FACE

- Create a curve in the flexible band for the nose bridge and expand the mask
- Secure the ties or elastics at the middle of the head and neck, or fasten the ear loops
- Fit it snugly to your face and below your chin – fit check the N95 and P2 mask
- Place the face shield or eyewear over the face and eyes and adjust to fit



3. PROTECT YOUR HANDS

- Extend the cuffs to cover the wrist of the gown
- Limit the surfaces touched
- Change gloves when torn or heavily contaminated
- Perform hand hygiene as required



FINALLY STOP and CHECK all PPE is in place



SEQUENCE FOR ADVANCED DOFFING PERSONAL PROTECTIVE EQUIPMENT (PPE):

Remember! The outside and front of your PPE, including the sleeves of the gown, are contaminated. If hands become visibly contaminated during PPE removal, wash hands before continuing.

REMOVE GOWN THEN GLOVES


Outside of gown and gloves are **CONTAMINATED!**

- 1a.**
- With gloved hands, grasp the front of the gown at the shoulders and pull the gown forward to snap the tabs at the back
 - Peel off the gown, rolling it inside out
- 1b.**
- Continue over the gloves so that they too are rolled inside out
 - Discard gown and gloves into a waste container
 - Perform hand hygiene using ABHR or soap and water



2. REMOVE PROTECTIVE EYEWEAR

- Grasp the ear or head piece and lift away from face
- Place reusable eyewear in the designated receptacle for cleaning or
- Discard disposable eyewear into waste container for disposal



3. REMOVE MASK OR PARTICULATE RESPIRATOR (N95 OR P2)

- Remove by touching ties or elastics only, starting with the bottom then top ties or bands
- Discard in designated waste container



4. PERFORM HAND HYGIENE

- Wash hands with soap and water or use an alcohol-based hand rub immediately after removing PPE



PROTECTIVE APPAREL



HALYARD* ISOLATION GOWNS

- For use when expected risk of exposure to fluid is between low and moderate
- Medium-weight tri-layer fabric
- Tape-tab neck closure
- Tie waist, elastic cuffs

| Code | Description | Packaging |
|-------------|---------------------------------|--|
| 69979 | Isolation Gown, Yellow, Large | Eaches/Bundle = 10 Bundles/Case = 10 Eaches/Case = 100 |
| 69988 | Isolation Gown, Yellow, X-Large | Eaches/Bundle = 10 Bundles/Case = 10 Eaches/Case = 100 |

PROTECTIVE APPAREL

HALYARD* IMPERVIOUS COMFORT GOWNS



69455

- Potential for moderate to heavy fluid contact in front
- Tri-layer spunbond with polyethylene laminate coating
- Breathable spunbond back
- Knit cuffs
- Velcro neck tabs
- Tie waist

| Code | Description | Packaging |
|-------|--|--|
| 69455 | Impervious Open Back Gown, Knit Cuff, Blue, Large | Eaches/Bundle = 30 Bundles/Case = 2 Eaches/Case = 60 |
| 5511 | Impervious Full Back Comfort Gown, Yellow, Regular | Eaches/Case = 50 |
| 5512 | Impervious Full Back Comfort Gown, Yellow, X-Large | Eaches/Case = 50 |



5511

IMPERVIOUS THUMBS-UP* FILM GOWNS

- Potential for soaking and heavy fluid contact in front
- Impervious polyethylene film
- Open back, thumbhooks



| Code | Description | Packaging |
|------|--|--|
| 7000 | Impervious Thumbs-Up* Film Gown with Thumbhooks, Blue, Regular | Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75 |
| 7001 | Impervious Thumbs-Up* Film Gown with Thumbhooks, Blue, X-Large | Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75 |
| 7002 | Impervious Thumbs-Up* Film Gown with Thumbhooks, Yellow, Regular | Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75 |
| 7003 | Impervious Thumbs-Up* Film Gown with Thumbhooks, Yellow, X-Large | Eaches/Dispenser = 15 Dispensers/Case = 5 Eaches/Case = 75 |

PROTECTIVE APPAREL



HALYARD* CHEMOTHERAPY PROCEDURE GOWN

- Potential for soaking and heavy fluid contact in front
- Tested against 52 chemotherapy drugs
- Meets Liquid Barrier Performance Standards ASTM 1670 and 1671
- Highest level of protection against flammability
- Polypropylene with polyethylene coating



| Code | Description | Packaging |
|-------|--|-------------------|
| 69606 | Procedure Gown for Chemotherapy Use, Over-the-head Open-back, Poly-coated, Universal | Eaches/Case = 100 |
| 37284 | Procedure Gown for Chemotherapy Use, Over-the-head Open-back, Poly-coated, XX-Large | Eaches/Case = 100 |



HALYARD* CHEMO360* GOWN

- Closed back chemotherapy gown
- Potential for soaking and heavy fluid contact in front
- Tested against 52 chemotherapy drugs
- Meets Liquid Barrier Performance Standards ASTM 1670 and 1671
- Highest level of protection against flammability
- Polypropylene with polyethylene coating



| Code | Description | Packaging |
|-------|---|-------------------|
| 47346 | CHEMO360* Gown, Closed-back, Poly-coated, Tape closure, Universal | Eaches/Case = 100 |
| 47347 | CHEMO360* Gown, Closed-back, Poly-coated, Tape closure, XX-Large | Eaches/Case = 100 |

PROTECTIVE APPAREL

SCRUB SUITS

HALYARD* SCRUB SHIRTS



- Potential for light fluid contact
- Protective tri-layer SMS fabric
- Scoop neck
- Reinforced pocket

| Code | Description | Packaging |
|-------------|---|------------------|
| 69701 | HALYARD* Scrub Shirt, Tri-Layer Fabric, Blue, Medium | Eaches/Case = 48 |
| 69702 | HALYARD* Scrub Shirt, Tri-Layer Fabric, Blue, Large | Eaches/Case = 48 |
| 69703 | HALYARD* Scrub Shirt, Tri-Layer Fabric, Blue, X-Large | Eaches/Case = 48 |

HALYARD* SCRUB PANTS

- Potential for light fluid contact
- Protective tri-layer SMS fabric
- Drawstring waist
- Side and rear pockets
- Reinforced in-seam

| Code | Description | Packaging |
|-------------|-------------------------------------|------------------|
| 69711 | HALYARD* Scrub Pants, Blue, Medium, | Eaches/Case = 48 |
| 69712 | HALYARD* Scrub Pants, Blue, Large | Eaches/Case = 48 |
| 69713 | HALYARD* Scrub Pants, Blue, X-Large | Eaches/Case = 48 |

PROTECTIVE APPAREL

HEADWEAR

HALYARD* BOUFFANT CAPS

- No potential fluid contact
- 1-Layer spunbond fabric
- Regular size



| Code | Description | Packaging |
|------|--------------------------------------|--|
| 3501 | Bouffant Cap, Spunbond Fabric, Blue | Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500 |
| 3502 | Bouffant Cap, Spunbond Fabric, White | Eaches/Dispenser = 100 Dispensers/Case = 5 Eaches/Case = 500 |

HALYARD* SURGICAL CAP

- For light fluid contact
- Tri-Layer SMS fabric on sides
- Spunbond fabric on top
- Universal size



| Code | Description | Packaging |
|-------|--------------------|--|
| 69520 | Surgical Cap, Blue | Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300 |

HALYARD* SURGICAL HOOD

- For light fluid contact
- Tri-layer SMS fabric on sides
- Spunbond fabric on top
- Tie neck
- Universal size



| Code | Description | Packaging |
|-------|--------------------------------|--|
| 69110 | Protective Surgical Hood, Blue | Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300 |

PROTECTIVE APPAREL

FOOTWEAR

X-TRA TRACTION* SHOE COVERS

- Blue
- Protective tri-layer SMS fabric
- Traction strips
- Natural rubber latex-free



| Code | Description | Packaging |
|-------|---------------------------------------|--|
| 69252 | X-TRA TRACTION* Shoe Cover, Universal | Eaches/Dispenser = 100 Dispensers/Case = 3 Eaches/Case = 300 |
| 69254 | X-TRA TRACTION* Shoe Cover, X-Large | Eaches/Dispenser = 80 Dispensers/Case = 3 Eaches/Case = 240 |

ANKLE TOPS

ANKLE GUARD* SHOE COVER

- Moderate fluid contamination
- Over the ankle coverage
- Protective tri-layer SMS fabric
- Adhesive traction strips
- Natural rubber latex-free



| Code | Description | Packaging |
|-------|--|---|
| 69353 | ANKLE GUARD* Shoe Cover, Tri-Layer Fabric, Adhesive Traction Strips, X-Large | Eaches/Dispenser = 90 Dispensers/Case = 3 Eaches/Case = 270 |

KNEE HIGHS

HI-GUARD* ULTRA FULL-COVERAGE BOOT

- Potential for heavy fluid contamination
- Knee high boot
- Plastic film coating over tri-layer fabric
- Two 1", 2.5cm adhesive foam traction strips
- Natural rubber latex-free



| Code | Description | Packaging |
|-------|---|---|
| 69572 | HI-GUARD* Full Coverage Boot, Universal | Eaches/Dispenser = 30 Dispensers/Case = 4 Eaches/Case = 120 |
| 69672 | HI-GUARD* Full Coverage Boot, X-Large | Eaches/Dispenser = 30 Dispensers/Case = 4 Eaches/Case = 120 |

NOTE: Please see package inserts for directions for use for each of these products. Always read the label and follow the directions for use.



Infection Prevention Product Solutions
Halyard* Education Foundation
Tools & Best Practice
Commitment to Excellence and Sustainability
Experienced Customer Support
Knowledgeable Sales Force
Innovation & Technology



O&M Halyard Australia Pty. Limited

AUS: 1800 664 227

www.halyardhealth.com.au

www.halyardhealth.com.au



Halyard makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2022 O&M Halyard, Inc.
All Rights Reserved. 2022:06:WFHY947v6