



# PROTECTIVE APPAREL

# PROTECTIVE APPAREL

## PROCEDURE AND ISOLATION GOWNS



### Procedure Gowns (Tested Against Chemotherapy Drugs)

- Tested for use with 52 chemotherapy drugs
- Polypropylene Gown with a Polyethylene Coating and Knit Cuffs
- Class 1 Flammability when tested as per NFPA Standard #702
- Potential for soaking and heavy fluid contact in front

| Code  | Description           | Items/Bundle | Bundles/Case |
|-------|-----------------------|--------------|--------------|
| 69606 | <b>PROCEDURE GOWN</b> | 10           | 10           |

Description: Knit Cuffs, Perforated Open Back, Blue, Large



### HALYARD\* Isolation Gowns

- Medium-Weight 3-Layer SMS Fabric
- For use when expected risk of exposure to fluid is between low and moderate

| Code  | Description                    | Items/Bundle | Bundles/Case |
|-------|--------------------------------|--------------|--------------|
| 69981 | <b>HALYARD* ISOLATION GOWN</b> | 10           | 10           |

Description: Elastic Cuff, Tie Waist, Tape-Tab Neck Closure, Blue, Large

|       |                                       |    |    |
|-------|---------------------------------------|----|----|
| 69987 | <b>TRI-LAYER AAMI2 ISOLATION GOWN</b> | 10 | 10 |
|-------|---------------------------------------|----|----|

Description: Elastic Cuff, Tie Waist, Tape-Tab Neck Closure, Blue, X-Large



### Impervious Thumbs-Up gown

- Impervious Polyethylene Film
- Potential for soaking and heavy fluid contact in front
- Packaged in a dispenser box

| Code | Description                 | Items/Bundle | Bundles/Case |
|------|-----------------------------|--------------|--------------|
| 7000 | <b>IMPERVIOUS FILM GOWN</b> | 15           | 5            |

Description: Thumbhooks, Open Back, Blue, Regular

# PROTECTIVE APPAREL

## HEADWEAR



### Surgical Caps

- Protective Surgical Cap
- 3-Layer Fabric

|                     |       |     |   |
|---------------------|-------|-----|---|
| <b>SURGICAL CAP</b> | 69540 | 100 | 3 |
|---------------------|-------|-----|---|

Protective 3-layer SMS fabric on side, Spunbond on top, Light fluid contact



### Surgical Hoods

- Protective 3-Layer SMS Fabric
- Tie Neck
- Universal Size

|                                 | <b>Code</b> | <b>Items/Dispenser</b> | <b>Dispensers/Case</b> |
|---------------------------------|-------------|------------------------|------------------------|
| <b>PROTECTIVE SURGICAL HOOD</b> | 69110       | 100                    | 3                      |

Description: Neck Tie, Blue, Universal

## FOOTWEAR



### ANKLE GUARD\* Shoe Cover

- Over the ankle coverage
- Protective 3-Layer SMS Fabric

|                                | <b>Code</b> | <b>Items/Dispenser</b> | <b>Dispensers/Case</b> |
|--------------------------------|-------------|------------------------|------------------------|
| <b>ANKLE GUARD* SHOE COVER</b> | 69353       | 90                     | 3                      |

Description: Adhesive Traction Strips, X-Large



### HI GUARD\* Ultra Full Coverage Boot

|   |       |    |   |
|---|-------|----|---|
| <b>HI GUARD* Ultra Full Coverage Boot</b> | 69572 | 30 | 4 |
|---|-------|----|---|

Description: Shoe and lower leg coverage  
Plastic film coating over SMS fabric; two 1" urethane foam traction strips  
Heavy fluid contact

# PROTECTIVE APPAREL


## SCRUB SUITS



### Scrub Shirts

Protective 3-Layer SMS Fabric  
Scoop Neck  
Reinforced Pocket  
Potential for light fluid contact

|                            | Code  | Items/Bundle | Bundles/Case |
|----------------------------|-------|--------------|--------------|
| <b>SCRUB SHIRT</b>         | 69701 | 1            | 48           |
| Description: Blue, Medium  |       |              |              |
| <b>SCRUB SHIRT</b>         | 69702 | 1            | 48           |
| Description: Blue, Large   |       |              |              |
| <b>SCRUB SHIRT</b>         | 69703 | 1            | 48           |
| Description: Blue, X-Large |       |              |              |



### Scrub Pants

Protective 3-Layer SMS Fabric  
Drawstring Waist  
Side and Rear Pockets  
Reinforced In-Seam  
Potential for light fluid contact

|                            |       |   |    |
|----------------------------|-------|---|----|
| <b>SCRUB PANTS</b>         | 69711 | 1 | 48 |
| Description: Blue, Medium  |       |   |    |
| <b>SCRUB PANTS</b>         | 69712 | 1 | 48 |
| Description: Blue, Large   |       |   |    |
| <b>SCRUB PANTS</b>         | 69713 | 1 | 48 |
| Description: Blue, X-Large |       |   |    |



**O&M Halyard India Pvt Ltd**  
World Trade Center, Tower-1 Unit 308  
Kharadi, Pune 411 014  
Maharashtra, India

For feedback/suggestions/complaints, contact Customer care officer  
Tel no: +91 20-46845500/5555  
Email: [customercare\\_IND@hyh.com](mailto:customercare_IND@hyh.com)

COPY-03816

O&M Halyard makes no representation that its products described in this catalogue will provide complete protection from exposure to contaminants, fluids or micro-organisms or will eliminate the risk of contracting infectious diseases.

\*Registered Trademark or Trademark of O&M Halyard Inc. or its affiliates.  
©2021. All rights reserved

[www.halyardhealth.in](http://www.halyardhealth.in)  
 O&M Halyard APAC