

# HALYARD\* PUREZERO\* NITRILE CLEANROOM GLOVES

HALYARD\* offers an accelerator-free<sup>1</sup> portfolio of three non-sterile and two sterile cleanroom gloves to address the needs of pharmaceutical, medical device, microelectronics and semi-conductor manufacturing industries.

## GLOVE SELECTION GUIDE

	Description	Designed for	Max Particle Count	ISO Class	Finish	Fingertips	Double Donning	Product Code	Fingertip Thickness	Case Count
NON - STERILE	HALYARD* PUREZERO* HG3 White Nitrile Gloves	Semiconductor, Pharmaceutical, Medical Device, Manufacturing,	<950	ISO Class 3 or higher and Grade A/B/C/D Cleanrooms	Tacky grip	Textured Fingertips	Recommended for Outer Glove	S: CLN3031SM M: CLN3031MD L: CLN3031LG	0.16mm (6 mil)	1000/cs (100X10) Ambidextrous
	HALYARD* PUREZERO* HG3 Light Blue Nitrile Gloves	Semiconductor, Pharmaceutical, Medical Device, Manufacturing,	<1200	ISO Class 3 or higher and Grade A/B/C/D Cleanrooms	Tacky grip	Textured Fingertips	Recommended for Outer Glove	S: CLN9031SM M: CLN9031MD L: CLN9031LG	0.10mm (4 mil)	1500/cs (250X6) Ambidextrous
	HALYARD* PUREZERO* HG3 SMOOTH White Nitrile Gloves	Semiconductor, Pharmaceutical, Medical Device, Manufacturing,	<950	ISO Class 3 or higher and Grade A/B/C/D Cleanrooms	Smooth grip	Textured Fingertips	Outer or Under	S: CLN3231SM M: CLN3231MD L: CLN3231LG	0.16mm (6 mil)	1000/cs (100X10) Ambidextrous
STERILE	HALYARD* PUREZERO* HG3 White Sterile Nitrile Gloves	Pharmaceutical, Biotechnology, Sterile Compounding, Aseptic Processing	<950	ISO Class 3 or higher and Grade A/B/C/D Cleanrooms	Smooth grip	Textured Fingertips and Palms	Outer or Under	6.5: CLN323265 7.0: CLN323270 7.5: CLN323275 8.0: CLN323280	0.16mm (6 mil)	200 pairs/cs (20 pairsX10) Hand Specific
	HALYARD* PUREZERO* HG3 Light Blue Sterile Nitrile Gloves	Pharmaceutical, Biotechnology, Sterile Compounding, Aseptic Processing	<1200	ISO Class 3 or higher and Grade A/B/C/D Cleanrooms	Smooth grip	Textured Fingertips and Palms	Outer or Under	6.5: CLN923265 7.0: CLN923270 7.5: CLN923275 8.0: CLN923280	0.10mm (4 mil)	300 pairs/cs (30 pairsX10) Hand Specific

Applies to all PUREZERO\* Gloves: Accelerator-Free<sup>1</sup> | 12" inch/30.5cm | AQL 1.0 | ISO 9001 | ISO13485 | **CE 2797** (PPE Cat. III) | ISO 374-5:2016 VIRUS | ISO 374-1:2016/Type C K-Low Chemical

For more information or samples, contact your distributor or visit: [www.halyardhealth.com/cleanroom](http://www.halyardhealth.com/cleanroom)



1 Not formulated with these commonly used vulcanizing chemicals: Sulfur, Thiurams, Thioxoles, Guanidines and Carbamates.  
 2 Tested per ASTM D6319, EN 455-2  
 \*Registered Trademark or Trademark of O&M Halyard or its affiliates. ©2020. All rights reserved.  
 H02757 H00078-12-01, 008271, COPY-03732